
ventilator

Release 0.0.0

jonny saunders et al

May 20, 2020

OVERVIEW

1	Vent control pseudocode	1
1.1	Pressure control parameters	1
1.2	Hardware	2
1.3	Pressure control loop	2
2	Hardware	5
3	common module	7
3.1	message	7
3.2	values	8
4	controller module	13
5	coordinator module	17
5.1	Submodules	17
5.2	coordinator	17
5.3	ipc	19
5.4	process_manager	19
6	gui	21
6.1	Program Diagram	21
6.2	Design Requirements	21
6.3	UI Notes & Todo	21
6.4	Jonny Questions	24
6.4.1	jonny todo	24
6.5	GUI Object Documentation	24
6.5.1	Control	24
6.5.2	Monitor	24
6.5.3	Plot	25
6.5.4	Status Bar	26
6.5.5	Components	28
6.5.5.1	GUI Stylesheets	28
6.5.5.2	GUI Alarm Manager	29
7	vent.io package	31
7.1	Subpackages	31
7.1.1	vent.io.devices package	31
7.1.1.1	Submodules	31
7.1.1.2	vent.io.devices.base module	31
7.1.1.3	vent.io.devices.pins module	40
7.1.1.4	vent.io.devices.sensors module	43

7.1.1.5	vent.io.devices.valves module	49
7.1.1.6	Module contents	55
7.2	Submodules	55
7.3	vent.io.hal module	55
7.4	Module contents	57
8	alarm	59
8.1	Main Alarm Module	59
8.2	Alarm Manager	60
8.3	Alarm	63
8.4	Condition	64
8.5	Alarm Rule	69
9	Requirements	71
10	Datasheets & Manuals	73
10.1	Manuals	73
10.2	Other Reference Material	73
11	Specs	75
12	Changelog	77
12.1	Version 0.0	77
12.1.1	v0.0.2 (April xxth, 2020)	77
12.1.2	v0.0.1 (April 12th, 2020)	77
12.1.3	v0.0.0 (April 12th, 2020)	77
13	Building the Docs	79
13.1	Local Build	79
14	h1 Heading 8-)	81
14.1	h2 Heading	81
14.1.1	h3 Heading	81
14.1.1.1	h4 Heading	81
14.2	Horizontal Rules	81
14.3	Emphasis	81
14.4	Blockquotes	81
14.5	Lists	82
14.6	Code	82
14.7	Links	83
14.8	Images	85
15	Indices and tables	87
	Python Module Index	89
	Index	91

VENT CONTROL PSEUDOCODE

Describes the procedure to operate low-cost ventilator under pressure control

1.1 Pressure control parameters

Set in GUI

- PIP: peak inhalation pressure (~20 cm H₂O)
- T_{insp}: inspiratory time to PEEP (~0.5 sec)
- I/E: inspiratory to expiratory time ratio
- bpm: breaths per minute (15 bpm -> 1/15 sec cycle time)
- PIP_{time}: Target time for PIP. While lungs expand, dP/dt should be PIP/PIP_{time}
- flow_{insp}: nominal flow rate during inspiration

Set by hardware

- FiO₂: fraction of inspired oxygen, set by blender
- max_flow: manual valve at output of blender
- PEEP: positive end-expiratory pressure, set by manual valve

Derived parameters

- cycle_time: 1/bpm
- t_{insp}: inspiratory time, controlled by cycle_time and I/E
- t_{exp}: expiratory time, controlled by cycle_time and I/E

Monitored variables

- Tidal volume: the volume of air entering the lung, derived from flow through t_{exp}
- PIP: peak inspiratory pressure, set by user in software
- Mean plateau pressure: derived from pressure sensor during inspiration cycle hold (no flow)
- PEEP: positive end-expiratory pressure, set by manual valve

Alarms

- Oxygen out of range
- High pressure (tube/airway occlusion)
- Low-pressure (disconnect)

- Temperature out of range
- Low voltage alarm (if using battery power)
- Tidal volume (expiratory) out of range

1.2 Hardware

Sensors

- O2 fraction, in inspiratory path
- Pressure, just before wye to endotracheal tube
- Flow, on expiratory path
- Temperature
- Humidity?

Actuators

- Inspiratory valve
 - Proportional or on/off
 - Must maintain low flow during expiratory cycle to maintain PEEP
- Expiratory valve
 - On/off in conjunction with PEEP valve probably OK

1.3 Pressure control loop

1. Begin inhalation
 - v1 Triggered by program every 1/bpm sec
 - v2 triggered by momentary drop in pressure when patient initiates inhalation (technically pressure-assisted control, PAC)
 - 1. ExpValve.close()
 - 2. InspValve.set(flow_insp)
2. While PSensor.read() < PIP
 - 1. Monitor $d(\text{PSensor.read()})/dt$
 - 2. Adjust flow rate for desired slope with controller
3. Cut flow and hold for t_{insp}
 - 1. InspValve.close()
 - 2. Monitor PSensor.read() and average across this time interval to report mean plateau pressure
4. Begin exhalation and hold for t_{exp}
 - 1. InspValve.set(PEEP_flow_rate) (alt: switch to parallel tube with continuous flow)
 - 2. ExpValve.open()
 - 3. integrate(FSensor.read()) for $t_{\text{exhalation}}$ to determine V_{tidal}

5. Repeat from step 1.

HARDWARE

COMMON MODULE

3.1 message

Classes

<code>ControlSetting(name, value, min_value, ...)</code>	TODO: if enum is hard to use, we may just use a predefined set, e.g.
<code>Error(errnum, err_str, timestamp)</code>	
<code>SensorValues([timestamp, loop_counter, ...])</code>	param timestamp from <code>time.time()</code> . must be passed explicitly or as an entry in <code>vals</code>

```
class vent.common.message.SensorValues (timestamp=None,          loop_counter=None,
                                         breath_count=None, vals=None, **kwargs)
```

Bases: `object`

Parameters

- **timestamp** (*float*) – from `time.time()`. must be passed explicitly or as an entry in `vals`
- **loop_counter** (*int*) – number of control_module loops. must be passed explicitly or as an entry in `vals`
- **breath_count** (*int*) – number of breaths taken. must be passed explicitly or as an entry in `vals`
- ****kwargs** – sensor readings, must be in `vent.values.SENSOR.keys`

Methods

<code>__init__([timestamp, loop_counter, ...])</code>	param timestamp from <code>time.time()</code> . must be passed explicitly or as an entry in <code>vals</code>
<code>to_dict()</code>	

Attributes

additional_values	Built-in immutable sequence.
-------------------	------------------------------

```
additional_values = ('timestamp', 'loop_counter', 'breath_count')
```

```
__init__(timestamp=None, loop_counter=None, breath_count=None, vals=None, **kwargs)
```

Parameters

- **timestamp** (*float*) – from time.time(). must be passed explicitly or as an entry in vals
- **loop_counter** (*int*) – number of control_module loops. must be passed explicitly or as an entry in vals
- **breath_count** (*int*) – number of breaths taken. must be passed explicitly or as an entry in vals
- ****kwargs** – sensor readings, must be in vent.values.SENSOR.keys

```
to_dict()
```

```
class vent.common.message.ControlSetting(name, value, min_value, max_value, timestamp)
```

Bases: `object`

TODO: if enum is hard to use, we may just use a predefined set, e.g. {'PIP', 'PEEP', ...} :param name: enum belong to ValueName :param value: :param min_value: :param max_value: :param timestamp: **Methods**

__init__(name, value, min_value, max_value, ...)	TODO: if enum is hard to use, we may just use a predefined set, e.g.
--	--

```
__init__(name, value, min_value, max_value, timestamp)
```

TODO: if enum is hard to use, we may just use a predefined set, e.g. {'PIP', 'PEEP', ...} :param name: enum belong to ValueName :param value: :param min_value: :param max_value: :param timestamp:

```
class vent.common.message.Error(errnum, err_str, timestamp)
```

Bases: `object`

3.2 values

Parameterization of variables and values

Data

CONTROL	Values to control but not monitor.
LIMITS	Values that are dependent on other values:
SENSOR	Values to monitor but not control.

Classes

Value(name, units, abs_range, safe_range, ...)	Definition of a value.
ValueName	An enumeration.

```
class vent.common.values.ValueName
```

Bases: `enum.Enum`

An enumeration. **Attributes**

BREATHS_PER_MINUTE	int([x]) -> integer
FIO2	int([x]) -> integer
HUMIDITY	int([x]) -> integer
IE_RATIO	int([x]) -> integer
INSPIRATION_TIME_SEC	int([x]) -> integer
PEEP	int([x]) -> integer
PEEP_TIME	int([x]) -> integer
PIP	int([x]) -> integer
PIP_TIME	int([x]) -> integer
PRESSURE	int([x]) -> integer
TEMP	int([x]) -> integer
VTE	int([x]) -> integer

```

PIP = 1
PIP_TIME = 2
PEEP = 3
PEEP_TIME = 4
BREATHS_PER_MINUTE = 5
INSPIRATION_TIME_SEC = 6
IE_RATIO = 7
FIO2 = 8
TEMP = 9
HUMIDITY = 10
VTE = 11
PRESSURE = 12

```

```

class vent.common.values.Value(name: str, units: str, abs_range: tuple, safe_range: tuple, decimals: int, control: bool, sensor: bool, default: (<class 'int'>, <class 'float'>) = None)

```

Bases: `object`

Definition of a value.

Used by the GUI and control module to set defaults.

Parameters

- **name** (*str*) – Human-readable name of the value
- **units** (*str*) – Human-readable description of units
- **abs_range** (*tuple*) – tuple of ints or floats setting the logical limit of the value, eg. a percent between 0 and 100, (0, 100)
- **safe_range** (*tuple*) – tuple of ints or floats setting the safe ranges of the value,

note:

this is not the same thing as the user-set alarm values, though the user-set alarm values are initialized as ``safe\_range``.

- **decimals** (*int*) – the number of decimals of precision used when displaying the value

Methods

<code>__init__(name, units, abs_range, safe_range, ...)</code>	Definition of a value.
<code>to_dict()</code>	

Attributes

<code>abs_range</code>
<code>control</code>
<code>decimals</code>
<code>default</code>
<code>name</code>
<code>safe_range</code>
<code>sensor</code>

`__init__` (*name*: str, *units*: str, *abs_range*: tuple, *safe_range*: tuple, *decimals*: int, *control*: bool, *sensor*: bool, *default*: (<class 'int'>, <class 'float'>) = None)
 Definition of a value.

Used by the GUI and control module to set defaults.

Parameters

- **name** (str) – Human-readable name of the value
- **units** (str) – Human-readable description of units
- **abs_range** (tuple) – tuple of ints or floats setting the logical limit of the value, eg. a percent between 0 and 100, (0, 100)
- **safe_range** (tuple) – tuple of ints or floats setting the safe ranges of the value,

note:

```
this is not the same thing as the user-set alarm values,
though the user-set alarm values are initialized as ``safe_
↪range``.
```

- **decimals** (int) – the number of decimals of precision used when displaying the value

property name

property abs_range

property safe_range

property decimals

property default

property control

property sensor

to_dict()

```
vent.common.values.SENSOR = OrderedDict([( <ValueName.FIO2: 8>, <vent.common.values.Value c
Values to monitor but not control.
```

Used to set alarms for out-of-bounds sensor values. These should be sent from the control module and not computed.:

```
{
    'name' (str): Human readable name,
    'units' (str): units string, (like degrees or %),
    'abs_range' (tuple): absolute possible range of values,
    'safe_range' (tuple): range outside of which a warning will be raised,
    'decimals' (int): The number of decimals of precision this number should be
    ↳displayed with
}
```

`vent.common.values.CONTROL = OrderedDict([(ValueName.PIP: 1), <vent.common.values.Value of ...>])`
 Values to control but not monitor.

Sent to control module to control operation of ventilator.:

```
{
    'name' (str): Human readable name,
    'units' (str): units string, (like degrees or %),
    'abs_range' (tuple): absolute possible range of values,
    'safe_range' (tuple): range outside of which a warning will be raised,
    'default' (int, float): the default value of the parameter,
    'decimals' (int): The number of decimals of precision this number should be
    ↳displayed with
}
```

`vent.common.values.LIMITS = {}`
 Values that are dependent on other values:

```
{
    "dependent_value": (
        ['value_1', 'value_2'],
        callable_returning_boolean
    )
}
```

Where the first argument in the tuple is a list of the values that will be given as argument to the `callable_returning_boolean` which will return whether (True) or not (False) a value is allowed.

CONTROLLER MODULE

Classes

Balloon_Simulator(leak)	This is a imple physics simulator for inflating a balloon.
ControlModuleBase()	This is an abstract class for controlling simulation and hardware.
ControlModuleDevice()	
ControlModuleSimulator()	

Functions

get_control_module([sim_mode])

class vent.controller.control_module.ControlModuleBase

Bases: `object`

This is an abstract class for controlling simulation and hardware.

1. All internal variables fall in three classes, denoted by the beginning of the variable:

- “COPY_varname”: These are copies (see 1.) that are regularly sync’ed with internal variables.
- “__varname”: These are variables only used in the ControlModuleBase-Class
- “_varname”: These are variables used in derived classes.

2. Internal variables should only to be accessed though the **set_** and **get_** functions. These functions act on COPIES of internal variables (“__” and “_”), that are sync’d every few iterations. How often this is done is adjusted by the variable self._NUMBER_CONTROLL_LOOPS_UNTIL_UPDATE. To avoid multiple threads manipulating the same variables at the same time, every manipulation of “**COPY_**” is surrounded by a thread lock.

Methods

_ControlModuleBase__analyze_last_waveform()	This goes through the last waveform, and updates VTE, PEEP, PIP, PIP_TIME, I_PHASE, FIRST_PEEP and BPM.
_ControlModuleBase__calculate_control_signal_in()	
_ControlModuleBase__get_PID_error(ytarget, ...)	
_ControlModuleBase__test_critical_level(x)	This tests whether a variable is within bounds.
_ControlModuleBase__update_alarms()	This goes through the values obtained from the last waveform, and updates alarms.

continues on next page

Table 3 – continued from previous page

<code>_PID_reset()</code>	resets the PID cycle to zero
<code>_PID_update(dt)</code>	This instantiates the control algorithms.
<code>_alarm_to_COPY()</code>	
<code>_controls_from_COPY()</code>	
<code>_get_control_signal_in()</code>	This is the PID controlled signal on the inspiratory side
<code>_get_control_signal_out()</code>	This is the control signal (open/close) on the expiratory side
<code>_initialize_set_to_COPY()</code>	
<code>_sensor_to_COPY()</code>	
<code>_start_mainloop()</code>	
<code>do_pid_control()</code>	
<code>do_state_control()</code>	
<code>get_active_alarms()</code>	
<code>get_control(control_setting_name)</code>	Gets values of the COPY of the control settings.
<code>get_logged_alarms()</code>	
<code>get_past_waveforms()</code>	
<code>get_sensors()</code>	
<code>get_target_waveform()</code>	
<code>is_running()</code>	
<code>set_control(control_setting)</code>	Updates the entry of COPY contained in the control settings
<code>start()</code>	
<code>stop()</code>	

Public Methods: `get_sensors()`: Returns a copy of the current sensor values. `get_alarms()`: Returns a List of all alarms, active and logged `get_active_alarms()`: Returns a Dictionary of all currently active alarms. `get_logged_alarms()`: Returns a List of logged alarms, up to maximum length of `self._RINGBUFFER_SIZE` `get_control(ControlSetting)`: Sets a controll-setting. Is updated at latest within `self._NUMBER_CONTROLL_LOOPS_UNTIL_UPDATE` `get_past_waveforms()`: Returns a List of waveforms of pressure and volume during at the last N breath cycles, $N < \text{self._RINGBUFFER_SIZE}$, AND clears this archive. `get_target_waveform()`: Returns a step-wise linear target waveform, as defined by the current settings. `start()`: Starts the main-loop of the controller `stop()`: Stops the main-loop of the controller

```

_initialize_set_to_COPY()
_alarm_to_COPY()
_sensor_to_COPY()
_controls_from_COPY()
get_sensors() → vent.common.message.SensorValues
get_active_alarms()
get_logged_alarms() → List[vent.alarm.alarm.Alarm]
set_control(control_setting: vent.common.message.ControlSetting)
    Updates the entry of COPY contained in the control settings
get_control(control_setting_name: vent.common.values.ValueName) →
    vent.common.message.ControlSetting
    Gets values of the COPY of the control settings.

```

`_get_control_signal_in()`

This is the PID controlled signal on the inspiratory side

`_get_control_signal_out()`

This is the control signal (open/close) on the expiratory side

`_PID_reset()`

resets the PID cycle to zero

`_PID_update(dt)`

This instantiates the control algorithms. During the breathing cycle, it goes through the four states:

- 1) Rise to PIP
- 2) Sustain PIP pressure
- 3) Quick fall to PEEP
- 4) Sustain PEEP pressure

Once the cycle is complete, it checks the cycle for any alarms, and starts a new one. A record of pressure/volume waveforms is kept in `self.__cycle_waveform_archive`

dt: Time since last update in seconds

`get_past_waveforms()`

`get_target_waveform()`

`_start_mainloop()`

`start()`

`stop()`

`is_running()`

`do_pid_control()`

`do_state_control()`

`_ControlModuleBase__analyze_last_waveform()`

This goes through the last waveform, and updates VTE, PEEP, PIP, PIP_TIME, I_PHASE, FIRST_PEEP and BPM.

`_ControlModuleBase__calculate_control_signal_in()`

`_ControlModuleBase__get_PID_error(ytarget, yis, dt)`

`_ControlModuleBase__test_critical_levels(min, max, value, name)`

This tests whether a variable is within bounds. If it is, and an alarm existed, then the “alarm_end_time” is set. If it is NOT, a new alarm is generated and appended to the alarm-list. Input:

min: minimum value (e.g. 2) max: maximum value (e.g. 5) value: test value (e.g. 3) name: parameter type (e.g. “PIP”, “PEEP” etc.)

`_ControlModuleBase__update_alarms()`

This goes through the values obtained from the last waveform, and updates alarms.

class `vent.controller.control_module.ControlModuleDevice`

Bases: `vent.controller.control_module.ControlModuleBase` **Methods**

`_sensor_to_COPY()`

`_start_mainloop()`

```
_sensor_to_COPY()
```

```
_start_mainloop()
```

```
class vent.controller.control_module.Balloon_Simulator(leak)
```

```
Bases: object
```

This is a imple physics simulator for inflating a balloon. For math, see https://en.wikipedia.org/wiki/Two-balloon_experiment **Methods**

<code>OUupdate(variable, dt, mu, sigma, tau)</code>	This is a simple function to produce an OU process.
<code>_reset()</code>	resets the ballon to standard parameters.
<code>get_pressure()</code>	
<code>get_volume()</code>	
<code>set_flow_in(Qin, dt)</code>	
<code>set_flow_out(Qout, dt)</code>	
<code>update(dt)</code>	

```
get_pressure()
```

```
get_volume()
```

```
set_flow_in(Qin, dt)
```

```
set_flow_out(Qout, dt)
```

```
update(dt)
```

```
OUupdate(variable, dt, mu, sigma, tau)
```

This is a simple function to produce an OU process. It is used as model for fluctuations in measurement variables. inputs: variable: float value at previous time step dt : timestep mu : mean sigma : noise amplitude tau : time scale returns: new_variable : value of “variable” at next time step

```
_reset()
```

resets the ballon to standard parameters.

```
class vent.controller.control_module.ControlModuleSimulator
```

```
Bases: vent.controller.control_module.ControlModuleBase Methods
```

<code>_ControlModuleSimulator__SimulatedPropValve(x, dt)</code>	This simulates the action of a proportional valve.
<code>_ControlModuleSimulator__SimulatedSolenoid(x)</code>	Depending on x, set flow to a binary value.
<code>_sensor_to_COPY()</code>	
<code>_start_mainloop()</code>	

```
_sensor_to_COPY()
```

```
_start_mainloop()
```

```
_ControlModuleSimulator__SimulatedPropValve(x, dt)
```

This simulates the action of a proportional valve. Flow-current-curve eye-balled from the datasheet of SMC PVQ31-5G-23-01N <https://www.ocpneumatics.com/content/pdfs/PVQ.pdf>

x: Input current [mA] dt: Time since last setting in seconds [for the LP filter]

```
_ControlModuleSimulator__SimulatedSolenoid(x)
```

Depending on x, set flow to a binary value. Here: flow is either 0 or 1l/sec

```
vent.controller.control_module.get_control_module(sim_mode=False)
```

COORDINATOR MODULE

5.1 Submodules

5.2 coordinator

Classes

`CoordinatorBase([sim_mode])`

`CoordinatorLocal([sim_mode])`

param `sim_mode`

`CoordinatorRemote([sim_mode])`

Functions

`get_coordinator([single_process, sim_mode])`

class `vent.coordinator.coordinator.CoordinatorBase` (*sim_mode=False*)

Bases: `object` **Methods**

`get_control(control_setting_name)`

`get_sensors()`

`is_running()`

`set_control(control_setting)`

`start()`

`stop()`

get_sensors () → `vent.common.message.SensorValues`

set_control (*control_setting: vent.common.message.ControlSetting*)

get_control (*control_setting_name: vent.common.values.ValueName*) →
`vent.common.message.ControlSetting`

start ()

is_running () → `bool`

stop ()

class `vent.coordinator.coordinator.CoordinatorLocal` (*sim_mode=False*)

Bases: `vent.coordinator.coordinator.CoordinatorBase`

Parameters `sim_mode` –**Methods**

<code>__init__([sim_mode])</code>	param <code>sim_mode</code>
<code>get_control(control_setting_name)</code>	
<code>get_sensors()</code>	
<code>is_running()</code>	Test whether the whole system is running
<code>set_control(control_setting)</code>	
<code>start()</code>	Start the coordinator.
<code>stop()</code>	Stop the coordinator.

`_is_running``.set()` when thread should stop**Type** `threading.Event``__init__ (sim_mode=False)`**Parameters `sim_mode` –****`_is_running`**`.set()` when thread should stop**Type** `threading.Event`**`get_sensors()`** → `vent.common.message.SensorValues`**`set_control`** (`control_setting: vent.common.message.ControlSetting`)**`get_control`** (`control_setting_name: vent.common.values.ValueName`) →
`vent.common.message.ControlSetting`**`start()`**

Start the coordinator. This does a soft start (not allocating a process).

`is_running()` → `bool`

Test whether the whole system is running

`stop()`

Stop the coordinator. This does a soft stop (not kill a process)

class `vent.coordinator.coordinator.CoordinatorRemote (sim_mode=False)`Bases: `vent.coordinator.coordinator.CoordinatorBase` **Methods**

<code>get_control(control_setting_name)</code>	
<code>get_sensors()</code>	
<code>is_running()</code>	Test whether the whole system is running
<code>set_control(control_setting)</code>	
<code>start()</code>	Start the coordinator.
<code>stop()</code>	Stop the coordinator.

`get_sensors()` → `vent.common.message.SensorValues`**`set_control`** (`control_setting: vent.common.message.ControlSetting`)**`get_control`** (`control_setting_name: vent.common.values.ValueName`) →
`vent.common.message.ControlSetting`

start()
Start the coordinator. This does a soft start (not allocating a process).

is_running() → bool
Test whether the whole system is running

stop()
Stop the coordinator. This does a soft stop (not kill a process)

`vent.coordinator.coordinator.get_coordinator` (*single_process=False*,
sim_mode=False) →
`vent.coordinator.coordinator.CoordinatorBase`

5.3 ipc

5.4 process_manager

Classes

`ProcessManager`(*sim_mode*[, *startCommandLine*,
...])

class `vent.coordinator.process_manager.ProcessManager` (*sim_mode*, *startCommandLine=None*, *maxHeartbeatInterval=None*)

Bases: `object` **Methods**

`heartbeat`(*timestamp*)
`restart_process`()
`start_process`()
`try_stop_process`()

start_process()
try_stop_process()
restart_process()
heartbeat (*timestamp*)

6.1 Program Diagram

6.2 Design Requirements

- Display Values
 - Value name, units, absolute range, recommended range, default range
 - VTE
 - FiO2
 - Humidity
 - Temperature
- Plots
- Controlled Values
 - PIP
 - PEEP
 - Inspiratory Time
- Value Dependencies

6.3 UI Notes & Todo

- Jonny add notes from helpful RT in discord!!!
- Top status Bar
 - Start/stop button
 - Status indicator - a clock that increments with heartbeats, or some other visual indicator that things are alright
 - Status bar - most recent alarm or notification w/ means of clearing
 - Override to give 100% oxygen and silence all alarms
- API
 - Two queues, input and output. Read from socket and put directly into queue.




Fig. 1: Schematic diagram of major UI components and signals

- Input, receive (timestamp, key, value) messages where key and value are names of variables and their values
- Output, send same format
- Menus
 - Trigger some testing/calibration routine
 - Log/alarm viewer
 - Wizard to set values?
 - save/load values
- Alarms
 - Multiple levels
 - Silenced/reset
 - Logging
 - Sounds?
- General
 - Reduce space given to waveforms
 - Clearer grouping & titling for display values & controls
 - Collapsible range setting
 - Ability to declare dependencies between values
 - * Limits - one value's range logically depends on another
 - * Derived - one value is computed from another/others
 - Monitored values should have defaults, warning range, and absolute range
 - Two classes of monitored values – ones with limits and ones without. There seem to be lots and lots of observed values, but only some need limits. might want to make larger drawer of smaller displayed values that don't need controls
 - Save/load parameters. Autosave, and autorestore if saved <5m ago, otherwise init from defaults.
 - Implement timed updates to plots to limit resource usage
 - Make class for setting values
 - Possible plots
 - * Pressure vs. flow
 - * flow vs volume
 - * volume vs time
- Performance
 - Cache drawText() calls in range selector by drawing to pixmap

6.4 Jonny Questions

- Which alarm sounds to use?
- Final final final breakdown on values and ranges plzzz

6.4.1 jonny todo

- use loop_counter to check on controller advancement
- implement single values list with properties ‘controllable’ vs not.

6.5 GUI Object Documentation

6.5.1 Control

6.5.2 Monitor

Classes

`Monitor(value[, update_period, enum_name])`

param value

```
class vent.gui.widgets.monitor.Monitor(value, update_period=0.5, enum_name=None)
```

Parameters

- **value** (Value) –
- **update_period** (*float*) – update period of monitor in s

Methods

`__init__(value[, update_period, enum_name])`

param value

`_limits_changed(val)`

`check_alarm([value])`

`init_ui()`

`set_alarm(alarm)`

Simple wrapper to set alarm state from a qt signal

`timed_update()`

`toggle_alarm()`

`toggle_control(state)`

`update_boxes(new_values)`

`update_limits(new_limits)`

`update_value(new_value)`

Attributes

alarm_state

```
alarm(*args, **kwargs) = <PySide2.QtCore.Signal object>
limits_changed(*args, **kwargs) = <PySide2.QtCore.Signal object>
__init__(value, update_period=0.5, enum_name=None)
```

Parameters

- **value** (Value) –
- **update_period** (*float*) – update period of monitor in s

```
init_ui()
toggle_control(state)
update_boxes(new_values)
update_value(new_value)
update_limits(new_limits)
timed_update()
_limits_changed(val)
property alarm_state
set_alarm(alarm)
    Simple wrapper to set alarm state from a qt signal
    Parameters alarm (bool) – Whether to set as alarm state or not
toggle_alarm()
check_alarm(value=None)
staticMetaObject = <PySide2.QtCore.QMetaObject object>
```

6.5.3 Plot

Data

PLOT_FREQ	Update frequency of Plot s in Hz
PLOT_TIMER	A QTimer that updates :class:`.TimedPlotCurveItem`s

Classes

Plot(name[, buffer_size, plot_duration, ...])	When initializing PlotWidget, <i>parent</i> and <i>background</i> are passed to GraphicsWidget.__init__() and all others are passed to PlotItem.__init__().
---	---

```
vent.gui.widgets.plot.PLOT_TIMER = None
    A QTimer that updates :class:`.TimedPlotCurveItem`s
vent.gui.widgets.plot.PLOT_FREQ = 5
    Update frequency of Plot s in Hz
```

class `vent.gui.widgets.plot.Plot` (*name*, *buffer_size*=4092, *plot_duration*=5, *abs_range*=None, *safe_range*=None, *color*=None, *units*=", ***kwargs*)

When initializing `PlotWidget`, *parent* and *background* are passed to `GraphicsWidget.__init__()` and all others are passed to `PlotItem.__init__()`. **Methods**

<code>_safe_limits_changed(val)</code>	
<code>set_duration(dur)</code>	
<code>set_safe_limits(limits)</code>	
<code>update_value(new_value)</code>	<i>new_value</i> : (timestamp from <code>time.time()</code> , <i>value</i>)

```
limits_changed(*args, **kwargs) = <PySide2.QtCore.Signal object>
set_duration(dur)
staticMetaObject = <PySide2.QtCore.QMetaObject object>
update_value(new_value)
    new_value: (timestamp from time.time(), value)
_safe_limits_changed(val)
set_safe_limits(limits)
```

6.5.4 Status Bar

Classes

<code>HeartBeat([update_interval, timeout_dur])</code>	param <i>update_interval</i> How often to do the heartbeat, in ms
--	--

<code>Message_Display()</code>
<code>Power_Button()</code>
<code>Status_Bar()</code>

- Start/stop button
-

class `vent.gui.widgets.status_bar.Status_Bar`

- Start/stop button
- **Status indicator - a clock that increments with heartbeats**, or some other visual indicator that things are alright
- Status bar - most recent alarm or notification w/ means of clearing
- Override to give 100% oxygen and silence all alarms

Methods

<code>init_ui()</code>

```
init_ui()
staticMetaObject = <PySide2.QtCore.QMetaObject object>
```

class `vent.gui.widgets.status_bar.Message_Display`

Attributes

alarm_level

Methods

clear_message()
draw_state([state])
init_ui()
make_icons()
update_message(alarm)

param alarm

message_cleared(*args, **kwargs) = <PySide2.QtCore.Signal object>

level_changed(*args, **kwargs) = <PySide2.QtCore.Signal object>

make_icons()

init_ui()

draw_state(state=None)

update_message(alarm)

Parameters alarm(Alarm) –

clear_message()

property alarm_level

staticMetaObject = <PySide2.QtCore.QMetaObject object>

class vent.gui.widgets.status_bar.HeartBeat(update_interval=100, timeout_dur=5000)

Parameters

- **update_interval** (*int*) – How often to do the heartbeat, in ms
- **timeout** (*int*) – how long to wait before hearing from control process

Methods

__init__([update_interval, timeout_dur])

param update_interval How often to
do the heartbeat, in ms

_heartbeat()

Called every (update_interval) milliseconds to set the
text of the timer.

beatheart(heartbeat_time)

check_timeout()

init_ui()

set_indicator([state])

start_timer([update_interval])

param update_interval How often (in
ms) the timer should be updated.

continues on next page

Table 11 – continued from previous page

<code>stop_timer()</code>	you can read the sign ya punk
<pre> timeout (*args, **kwargs) = <PySide2.QtCore.Signal object> heartbeat (*args, **kwargs) = <PySide2.QtCore.Signal object> __init__ (update_interval=100, timeout_dur=5000) </pre>	
Parameters	
• update_interval (<i>int</i>) – How often to do the heartbeat, in ms • timeout (<i>int</i>) – how long to wait before hearing from control process	
<pre> init_ui() check_timeout() set_indicator (state=None) start_timer (update_interval=None) </pre>	
Parameters update_interval (<i>float</i>) – How often (in ms) the timer should be updated.	
<pre> stop_timer() you can read the sign ya punk beatheart (heartbeat_time) _heartbeat() Called every (update_interval) milliseconds to set the text of the timer. staticMetaObject = <PySide2.QtCore.QMetaObject object> </pre>	
class <code>vent.gui.widgets.status_bar.Power_Button</code>	
Methods	
<code>init_ui()</code>	
<pre> init_ui() staticMetaObject = <PySide2.QtCore.QMetaObject object> </pre>	

6.5.5 Components

6.5.5.1 GUI Stylesheets

Data

<code>MONITOR_UPDATE_INTERVAL</code>	inter-update interval (seconds) for Monitor
--------------------------------------	---

Functions

<code>set_dark_palette(app)</code>	Apply Dark Theme to the Qt application instance.
------------------------------------	--

```

vent.gui.styles.MONITOR_UPDATE_INTERVAL = 0.5
    inter-update interval (seconds) for Monitor

```

Type (float)

`vent.gui.styles.set_dark_palette(app)`
Apply Dark Theme to the Qt application instance.

borrowed from <https://github.com/gmarull/qtmodern/blob/master/qtmodern/styles.py>

Args: `app` (QApplication): QApplication instance.

6.5.5.2 GUI Alarm Manager

Classes

`AlarmManager()`

Functions

`get_alarm_manager()`

`vent.gui.alarm_manager.get_alarm_manager()`

class `vent.gui.alarm_manager.AlarmManager`

Methods

<code>monitor_alarm(alarm)</code>	Parse a tentative alarm from a monitor – we should have already gotten an alarm from the controller, so this largely serves as a double check.
<code>parse_message(alarm)</code>	If an alarm doesn't have a message attr, make one for it.
<code>update_alarms(alarms)</code>	

new_alarm (*args, **kwargs) = <PySide2.QtCore.Signal object>

update_alarms (alarms)

monitor_alarm (alarm)

Parse a tentative alarm from a monitor – we should have already gotten an alarm from the controller, so this largely serves as a double check.

Doesn't use the `Alarm` class because creating a new alarm increments the counter.

Parameters `alarm` (tuple) – (monitor_name, monitor_value, timestamp)

parse_message (alarm)

If an alarm doesn't have a message attr, make one for it.

staticMetaObject = <PySide2.QtCore.QMetaObject object>

VENT.IO PACKAGE

7.1 Subpackages

7.1.1 vent.io.devices package

7.1.1.1 Submodules

7.1.1.2 vent.io.devices.base module

Base classes & functions used throughout vent.io.devices

Classes

<i>ADS1015</i> ([address, i2c_bus, pig])	ADS1015 16 bit, 4 Channel Analog to Digital Converter.
<i>ADS1115</i> ([address, i2c_bus, pig])	ADS1115 16 bit, 4 Channel Analog to Digital Converter.
<i>I2CDevice</i> (i2c_address, i2c_bus[, pig])	A class wrapper for pigpio I2C handles.
<i>IODeviceBase</i> (pig)	Abstract base Class for pigpio handles (or whatever other GPIO library
<i>SPIDevice</i> (channel, baudrate[, pig])	A class wrapper for pigpio SPI handles.

Functions

<i>be16_to_native</i> (data[, signed])	Unpacks a bytes-like object respecting big-endianness of outside world and returns an int according to signed.
<i>native16_to_be</i> (word[, signed])	Packs an int into bytes after swapping endianness.

```
class vent.io.devices.base.IODeviceBase (pig:    vent.io.devices.base.PigpioConnection    =  
                                           None)
```

Bases: `object`

Abstract base Class for pigpio handles (or whatever other GPIO library we end up using)

Note: pigpio commands return -144 if an error is encountered while attempting to communicate with the demon. TODO would be to recognize when that occurs and handle it gracefully, i.e. kill the daemon, restart it, and reopen the python interface(s)

Initializes the pigpio python bindings object if necessary, and checks that it is actually running.

Parameters **pig** (*PigpioConnection*) – pigpiod connection to use; if not specified, a new one is established

Methods

<code>__init__(pig)</code>	Initializes the pigpio python bindings object if necessary, and checks that it is actually running.
<code>_close()</code>	Closes an I2C/SPI (or potentially Serial) connection

Attributes

<code>handle</code>	Pigpiod handle associated with device (only for i2c/spi)
<code>pig</code>	The pigpio python bindings object
<code>pigpiod_ok</code>	Returns True if pigpiod is running and False if not

`__init__` (*pig: vent.io.devices.base.PigpioConnection = None*)

Initializes the pigpio python bindings object if necessary, and checks that it is actually running.

Parameters **pig** (*PigpioConnection*) – pigpiod connection to use; if not specified, a new one is established

property pig

The pigpio python bindings object

property handle

Pigpiod handle associated with device (only for i2c/spi)

property pigpiod_ok

Returns True if pigpiod is running and False if not

`_close()`

Closes an I2C/SPI (or potentially Serial) connection

class `vent.io.devices.base.I2CDevice` (*i2c_address, i2c_bus, pig=None*)

Bases: `vent.io.devices.base.IODeviceBase`

A class wrapper for pigpio I2C handles. Defines several methods used for reading from and writing to device registers. Defines helper classes Register and ValueField for handling the manipulation of arbitrary registers.

Note: The Raspberry Pi uses LE byte-ordering, while the outside world tends to use BE (at least, the sensors in use so far all do). Thus, bytes need to be swapped from native (LE) ordering to BE prior to being written to an i2c device, and bytes recieved need to be swapped from BE into native (LE). All methods except `read_device` and `write_device` perform this automatically. The methods `read_device` and `write_device` do NOT byteswap and return bytearrays rather than the unsigned 16-bit int used by the other read/write methods.

Initializes pigpio bindings and opens i2c connection.

Parameters

- **i2c_address** (*int*) – I2C address of the device. (e.g., *i2c_address=0x50*)
- **i2c_bus** (*int*) – The I2C bus to use. Should probably be set to 1 on Raspberry Pi.
- **pig** (*PigpioConnection*) – pigpiod connection to use; if not specified, a new one is established

Classes

<code>Register</code> (fields, values)	Describes a writable configuration register.
--	--

Methods

<code>__init__(i2c_address, i2c_bus[, pig])</code>	Initializes pigpio bindings and opens i2c connection.
<code>_close()</code>	Extends superclass method.
<code>_open(i2c_bus, i2c_address)</code>	Opens i2c connection given i2c bus and address.
<code>read_device([count])</code>	Read a specified number of bytes directly from the the device without specifying or changing the register.
<code>read_register(register[, signed])</code>	Read 2 bytes from the specified register and byteswap the result.
<code>write_device(word[, signed])</code>	Write 2 bytes to the device without specifying register.
<code>write_register(register, word[, signed])</code>	Write 2 bytes to the specified register.

`__init__ (i2c_address, i2c_bus, pig=None)`
Initializes pigpio bindings and opens i2c connection.

Parameters

- **i2c_address** (*int*) – I2C address of the device. (e.g., *i2c_address=0x50*)
- **i2c_bus** (*int*) – The I2C bus to use. Should probably be set to 1 on Raspberry Pi.
- **pig** (*PigpioConnection*) – pigpiod connection to use; if not specified, a new one is established

`_open (i2c_bus, i2c_address)`
Opens i2c connection given i2c bus and address.

`_close ()`
Extends superclass method. Checks that pigpiod is connected and if a handle has been set - if so, closes an i2c connection.

`read_device (count=2) → tuple`
Read a specified number of bytes directly from the the device without specifying or changing the register. Does NOT perform LE/BE conversion.

Parameters **count** (*int*) – The number of bytes to read from the device.

Returns a tuple of the number of bytes read and a bytearray containing the bytes. If there was an error the number of bytes read will be less than zero (and will contain the error code).

Return type *tuple*

`write_device (word, signed=False)`
Write 2 bytes to the device without specifying register. DOES perform LE/BE conversion.

Parameters

- **word** (*int*) – The integer representation of the data to write.
- **signed** (*bool*) – Whether or not *word* is signed.

`read_register (register, signed=False) → int`
Read 2 bytes from the specified register and byteswap the result.

Parameters

- **register** (*int*) – The index of the register to read.
- **signed** (*bool*) – Whether or not the data to read is expected to be signed.

Returns integer representation of 16 bit register contents.

Return type *int*

write_register (*register*, *word*, *signed=False*)

Write 2 bytes to the specified register. Byteswaps.

Parameters

- **register** (*int*) – The index of the register to write to
- **word** (*int*) – The unsigned 16 bit integer to write to the register (must be consistent with ‘signed’)
- **signed** (*bool*) – Whether or not ‘word’ is signed

class Register (*fields*, *values*)

Bases: `object`

Describes a writable configuration register. Has dynamically defined attributes corresponding to the fields described by the passed arguments. Takes as arguments two tuples of equal length, the first of which names each field and the second being a tuple of tuples containing the (human readable) possible settings & values for each field.

Note: The initializer reverses the fields & their values because a human reads the register, as drawn in the datasheet, from left to right - however, the fields furthest to the left are the most significant bits of the register.

Initializer which loads (dynamically defined) attributes from tuples.

Parameters

- **fields** (*tuple*) – A tuple containing the names of the register’s value fields
- **values** (*tuple*) – A tuple of tuples containing the possible values for each value field. Length must match the length of fields. If there are redundant values for a field specified in the datasheet, be sure to include them. (e.g., a field takes values *A: 0b00*, *B: 0b01*, and *C: 0b10*; but the value for *0b11* is either not specified by the datasheet or is listed redundantly as *C: 0b11* -> *values* should list both the 3rd and 4th possible values as ‘C’ like so: (‘A’, ‘B’, ‘C’, ‘C’)

__init__ (*fields*, *values*)

Initializer which loads (dynamically defined) attributes from tuples.

Parameters

- **fields** (*tuple*) – A tuple containing the names of the register’s value fields
- **values** (*tuple*) – A tuple of tuples containing the possible values for each value field. Length must match the length of fields. If there are redundant values for a field specified in the datasheet, be sure to include them. (e.g., a field takes values *A: 0b00*, *B: 0b01*, and *C: 0b10*; but the value for *0b11* is either not specified by the datasheet or is listed redundantly as *C: 0b11* -> *values* should list both the 3rd and 4th possible values as ‘C’ like so: (‘A’, ‘B’, ‘C’, ‘C’)

unpack (*cfg*) → `collections.OrderedDict`

Given the contents of a register in integer form, returns a dict of fields and their current settings.

Parameters **cfg** (*int*) – An integer representing a possible configuration value for the register

pack (*cfg*, ***kwargs*) → *int*

Given an initial integer representation of a register and an arbitrary number of field=value settings, returns an integer representation of the register incorporating the new settings.

Parameters

- **cfg** (*int*) – An integer representing a possible configuration value for the register
- ****kwargs** – The register fields & values to patch into *cfg*. Takes keyword arguments of the form: *field=value*

class ValueField(*offset, mask, values*)

Bases: `object`

Describes a configurable value field in a writable register.

Instantiates a value field of a register given the bit offset, mask, and list of possible values.

Parameters

- **offset** (*int*) – The offset bits of the value field in the register, i.e. the distance from LSB
- **mask** (*int*) – integer representation of the value field mask (w/o offset)
- **values** (*OrderedDict*) – The possible values that the field can take.

___init___(*offset, mask, values*)

Instantiates a value field of a register given the bit offset, mask, and list of possible values.

Parameters

- **offset** (*int*) – The offset bits of the value field in the register, i.e. the distance from LSB
- **mask** (*int*) – integer representation of the value field mask (w/o offset)
- **values** (*OrderedDict*) – The possible values that the field can take.

unpack (*cfg*)

Extracts the ValueField's setting from *cfg* & returns the result in a human readable form.

Parameters **cfg** (*int*) – An integer representing a possible configuration value for the register

extract (*cfg*) → *int*

Extracts setting from passed 16-bit config & returns integer representation.

Parameters **cfg** (*int*) – An integer representing a possible configuration value for the register

pack (*value*) → *int*

Takes a human-readable ValueField setting and returns the corresponding bit-shifted integer.

Parameters **value** (*int*) – The integer representation of *value* bit-shifted by the ValueField's offset

Returns The integer representation of the ValueField setting according to *value*

Return type *int*

insert (*cfg, value*) → *int*

Validates and performs bitwise replacement with the human-readable ValueField setting and integer representation of the register configuration.

Parameters

- **cfg** (*int*) – An integer representing a possible configuration value for the register
- **value** (*object*) – The human readable representation of the desired ValueField setting. Must match a value in ValueField._values; if not, throws a ValueError

Returns

The integer representation of the Register's configuration with the value of ValueField patched according the *value*

Return type *int*

class `vent.io.devices.base.SPIDevice` (*channel, baudrate, pig=None*)

Bases: `vent.io.devices.base.IODeviceBase`

A class wrapper for pigpio SPI handles. Not really implemented.

Instantiates an SPIDevice on SPI *channel* with *baudrate* and, optionally, *pigpio.pi = pig*.

Parameters

- **channel** (*int*) – The SPI channel
- **baudrate** (*int*) – SPI baudrate

- **pig** (*PigpioConnection*) – pigpiod connection to use; if not specified, a new one is established

Methods

<code>__init__(channel, baudrate[, pig])</code>	Instantiates an SPIDevice on SPI <i>channel</i> with <i>baudrate</i> and, optionally, <i>pigpio.pi = pig</i> .
<code>__close()</code>	Extends superclass method.
<code>__open(channel, baudrate)</code>	Opens an SPI connection and sets the pigpiod handle.

`__init__ (channel, baudrate, pig=None)`

Instantiates an SPIDevice on SPI *channel* with *baudrate* and, optionally, *pigpio.pi = pig*.

Parameters

- **channel** (*int*) – The SPI channel
- **baudrate** (*int*) – SPI baudrate
- **pig** (*PigpioConnection*) – pigpiod connection to use; if not specified, a new one is established

`__open (channel, baudrate)`

Opens an SPI connection and sets the pigpiod handle.

Parameters

- **channel** (*int*) – The SPI channel
- **baudrate** (*int*) – SPI baudrate

`__close ()`

Extends superclass method. Checks that pigpiod is connected and if a handle has been set - if so, closes an SPI connection.

class `vent.io.devices.base.ADS1115 (address=72, i2c_bus=1, pig=None)`

Bases: `vent.io.devices.base.I2CDevice`

ADS1115 16 bit, 4 Channel Analog to Digital Converter. Datasheet:

Attributes

<code>USER_CONFIGURABLE_FIELDS</code>	Note: The Conversion Register is read-only and contains a 16bit representation of the requested value (provided the conversion is ready).
<code>__CONFIG_FIELDS</code>	Built-in immutable sequence.
<code>__CONFIG_VALUES</code>	Built-in immutable sequence.
<code>__DEFAULT_ADDRESS</code>	<code>int([x]) -> integer</code>
<code>__DEFAULT_VALUES</code>	<code>dict() -> new empty dictionary</code>
<code>__POINTER_FIELDS</code>	Built-in immutable sequence.
<code>__POINTER_VALUES</code>	Config Register (R/W)
<code>__TIMEOUT</code>	Address Pointer Register (write-only)
<code>cfg</code>	Returns the contents (as a 16-bit unsigned integer) of the configuration that will be written to the config register when <code>read_conversion()</code> is next called.
<code>config</code>	Returns the Register object of the config register.

Methods

<code>__init__([address, i2c_bus, pig])</code>	Initializes registers: Pointer register is write only, config is R/W.
<code>_read_conversion(**kwargs)</code>	Backend for read_conversion.
<code>_read_last_cfg()</code>	Reads the config register and returns the contents as a 16-bit unsigned integer; updates internal record <code>_last_cfg</code> .
<code>_ready()</code>	Return status of ADC conversion; True indicates the conversion is complete and the results ready to be read.
<code>print_config()</code>	Returns the human-readable configuration for the next read.
<code>read_conversion(**kwargs)</code>	Returns a voltage (expressed as a float) corresponding to a channel on the ADC.

<http://www.ti.com/lit/ds/symlink/ads1114.pdf?ts=1587872241912>

Default Values: Default configuration for vent: 0xC3E3 Default configuration on power-up: 0x8583

Initializes registers: Pointer register is write only, config is R/W. Sets initial value of `_last_cfg` to what is actually on the ADS. Packs default settings into `_cfg`, but does not actually write to ADC - that occurs when `read_conversion()` is called.

Parameters

- **address** (*int*) – I2C address of the device. (e.g., `i2c_address=0x48`)
- **i2c_bus** (*int*) – The I2C bus to use. Should probably be set to 1 on Raspberry Pi.
- **pig** (*PigpioConnection*) – pigpiod connection to use; if not specified, a new one is established

`_DEFAULT_ADDRESS = 72`

`_DEFAULT_VALUES = {'DR': 860, 'MODE': 'SINGLE', 'MUX': 0, 'PGA': 4.096}`

`_TIMEOUT = 1`

Address Pointer Register (write-only)

`_POINTER_FIELDS = ('P',)`

`_POINTER_VALUES = (('CONVERSION', 'CONFIG', 'LO_THRESH', 'HIGH_THRESH'),)`

Config Register (R/W)

`_CONFIG_FIELDS = ('OS', 'MUX', 'PGA', 'MODE', 'DR', 'COMP_MODE', 'COMP_POL', 'COMP_LAT')`

`_CONFIG_VALUES = (('NO_EFFECT', 'START_CONVERSION'), ((0, 1), (0, 3), (1, 3), (2, 3),`

`USER_CONFIGURABLE_FIELDS = ('MUX', 'PGA', 'MODE', 'DR')`

Note: The Conversion Register is read-only and contains a 16bit representation of the requested value (provided the conversion is ready).

The Lo-thresh & Hi-thresh Registers are not Utilized here. However, their function and usage are described in the datasheet. Should you want to extend the functionality implemented here.

`__init__(address=72, i2c_bus=1, pig=None)`

Initializes registers: Pointer register is write only, config is R/W. Sets initial value of `_last_cfg` to what is actually on the ADS. Packs default settings into `_cfg`, but does not actually write to ADC - that occurs when `read_conversion()` is called.

Parameters

- **address** (*int*) – I2C address of the device. (e.g., *i2c_address=0x48*)
- **i2c_bus** (*int*) – The I2C bus to use. Should probably be set to 1 on Raspberry Pi.
- **pig** (*PigpioConnection*) – pigpiod connection to use; if not specified, a new one is established

read_conversion (***kwargs*) → *float*

Returns a voltage (expressed as a float) corresponding to a channel on the ADC. The channel to read from, along with the gain, mode, and sample rate of the conversion may be specified as optional parameters. If `read_conversion()` is called with no parameters, the resulting voltage corresponds to the channel last read from and the same conversion settings.

Parameters

- **MUX** – The pin to read from in single channel mode: e.g., *0, 1, 2, 3* or, a tuple of pins over which to make a differential reading. e.g., *(0, 1), (0, 3), (1, 3), (2, 3)*
- **PGA** – The full scale voltage (FSV) corresponding to a programmable gain setting. e.g., *(6.144, 4.096, 2.048, 1.024, 0.512, 0.256, 0.256, 0.256)*
- **MODE** – Whether to set the ADC to continuous conversion mode, or operate in single-shot mode. e.g., *'CONTINUOUS', 'SINGLE'*
- **DR** – The data rate to make the conversion at; units: samples per second. e.g., *8, 16, 32, 64, 128, 250, 475, 860*

print_config () → *collections.OrderedDict*

Returns the human-readable configuration for the next read.

Returns an ordered dictionary of the form {field: value}, ordered from MSB -> LSB

Return type *OrderedDict*

property config

Returns the Register object of the config register.

Returns The Register object initialized for the ADS1115.

Return type *vent.io.devices.I2CDevice.Register*

property cfg

Returns the contents (as a 16-bit unsigned integer) of the configuration that will be written to the config register when `read_conversion()` is next called.

_read_conversion (***kwargs*) → *int*

Backend for `read_conversion`. Returns the contents of the 16-bit conversion register as an unsigned integer.

If no parameters are passed, one of two things can happen:

- 1) If the ADC is in single-shot (mode='SINGLE') conversion mode, `_last_cfg` is written to the config register; once the ADC indicates it is ready, the contents of the conversion register are read and the result is returned.
- 2) If the ADC is in CONTINUOUS mode, the contents of the conversion register are read immediately and returned.

If any of channel, gain, mode, or data_rate are specified as parameters, a new `_cfg` is packed and written to the config register; once the ADC indicates it is ready, the contents of the conversion register are read and the result is returned.

Note: In continuous mode, data can be read from the conversion register of the ADS1115 at any time and always reflects the most recently completed conversion. So says the datasheet.

Parameters ****kwargs** – see documentation of `vent.io.devices.ADS1115.read_conversion`

`_read_last_cfg()` → `int`

Reads the config register and returns the contents as a 16-bit unsigned integer; updates internal record `_last_cfg`.

`_ready()` → `bool`

Return status of ADC conversion; True indicates the conversion is complete and the results ready to be read.

class `vent.io.devices.base.ADS1015` (`address=72`, `i2c_bus=1`, `pig=None`)

Bases: `vent.io.devices.base.ADS1115`

ADS1015 16 bit, 4 Channel Analog to Digital Converter. Datasheet:

Attributes

<code>USER_CONFIGURABLE_FIELDS</code>	Built-in immutable sequence.
<code>_CONFIG_FIELDS</code>	Built-in immutable sequence.
<code>_CONFIG_VALUES</code>	Built-in immutable sequence.
<code>_DEFAULT_ADDRESS</code>	<code>int([x]) -> integer</code>
<code>_DEFAULT_VALUES</code>	Address Pointer Register (write-only)
<code>_POINTER_FIELDS</code>	Built-in immutable sequence.
<code>_POINTER_VALUES</code>	Config Register (R/W)

Methods

<code>__init__</code> (<code>[address, i2c_bus, pig]</code>)	See: <code>vent.io.devices.ADS1115.__init__</code>
--	--

<http://www.ti.com/lit/ds/symlink/ads1015.pdf?ts=1589228228921>

Basically the same device as the ADS1115, except has 12 bit resolution instead of 16, and has different (faster) data rates. The difference in data rates is handled by overloading `_CONFIG_VALUES`. The difference in resolution is irrelevant for implementation.

See: `vent.io.devices.ADS1115.__init__`

`_DEFAULT_ADDRESS = 72`

`_DEFAULT_VALUES = {'DR': 3300, 'MODE': 'SINGLE', 'MUX': 0, 'PGA': 4.096}`

Address Pointer Register (write-only)

`_POINTER_FIELDS = ('P',)`

`_POINTER_VALUES = (('CONVERSION', 'CONFIG', 'LO_THRESH', 'HIGH_THRESH'),)`

Config Register (R/W)

`_CONFIG_FIELDS = ('OS', 'MUX', 'PGA', 'MODE', 'DR', 'COMP_MODE', 'COMP_POL', 'COMP_LAT')`

`_CONFIG_VALUES = (('NO_EFFECT', 'START_CONVERSION'), ((0, 1), (0, 3), (1, 3), (2, 3),`

`USER_CONFIGURABLE_FIELDS = ('MUX', 'PGA', 'MODE', 'DR')`

`__init__` (`address=72`, `i2c_bus=1`, `pig=None`)

See: `vent.io.devices.ADS1115.__init__`

`vent.io.devices.base.be16_to_native` (`data`, `signed=False`) → `int`

Unpacks a bytes-like object respecting big-endianness of outside world and returns an int according to signed.

Parameters

- **data** – bytes-like object. The data to be unpacked & converted
- **signed** (*bool*) – Whether or not *data* is signed

`vent.io.devices.base.native16_to_be(word, signed=False) → bytes`

Packs an int into bytes after swapping endianness.

Parameters

- **signed** (*bool*) – Whether or not *data* is signed
- **word** (*int*) – The integer representation to converted and packed into bytes

7.1.1.3 vent.io.devices.pins module

Classes

<code>PWMOutput(pin[, initial_duty, frequency, pig])</code>	A pin configured to output a PWM signal.
<code>Pin(pin[, pig])</code>	Base Class wrapping pigpio methods for interacting with GPIO pins on the raspberry pi.

class `vent.io.devices.pins.Pin` (*pin, pig=None*)

Bases: `vent.io.devices.base.IODeviceBase`

Base Class wrapping pigpio methods for interacting with GPIO pins on the raspberry pi. Subclasses include InputPin, OutputPin; along with any specialized pins or specific devices defined in `vent.io.actuators` & `vent.io.sensors` (note: actuators and sensors do not need to be tied to a GPIO pin and may instead be interfaced through an ADC or I2C).

This is an abstract base class. The subclasses InputPin and OutputPin extend Pin into a usable form.

Inherits attributes and methods from IODeviceBase.

Parameters

- **pin** (*int*) – The number of the pin to use
- **pig** (*PigpioConnection*) – pigpiod connection to use; if not specified, a new one is established

Attributes

<code>__PIGPIO_MODES</code>	dict() -> new empty dictionary
<code>mode</code>	The currently active pigpio mode of the pin.

Methods

<code>__init__(pin[, pig])</code>	Inherits attributes and methods from IODeviceBase.
<code>read()</code>	Returns the value of the pin: usually 0 or 1 but can be overridden by subclass.
<code>toggle()</code>	If pin is on, turn it off.
<code>write(value)</code>	Sets the value of the Pin.

```
__PIGPIO_MODES = {'ALT0': 4, 'ALT1': 5, 'ALT2': 6, 'ALT3': 7, 'ALT4': 3, 'ALT5': 2}
```

```
__init__(pin, pig=None)
```

Inherits attributes and methods from IODeviceBase.

Parameters

- **pin** (*int*) – The number of the pin to use
- **pig** (*PigpioConnection*) – pigpiod connection to use; if not specified, a new one is established

property mode

The currently active pigpio mode of the pin.

toggle()

If pin is on, turn it off. If it's off, turn it on. Do not raise a warning when pin is read in this way.

read() → *int*

Returns the value of the pin: usually 0 or 1 but can be overridden by subclass.

write(value)

Sets the value of the Pin. Usually 0 or 1 but behavior differs for some subclasses.

Parameters value – The value to write to the pin. Can be either *1* to turn on the pin or *0* to turn it off.

class `vent.io.devices.pins.PWMOutput` (*pin*, *initial_duty=0*, *frequency=None*, *pig=None*)

Bases: `vent.io.devices.pins.Pin`

A pin configured to output a PWM signal. Can be configured to use either a hardware-generated or software-generated signal. Overrides parent methods `read()` and `write()`.

Inherits attributes from parent `Pin`, then sets PWM frequency & initial duty (use defaults if `None`)

Parameters

- **pin** (*int*) – The number of the pin to use. Hardware PWM pins are 12, 13, 18, and 19.
- **initial_duty** (*float*) – The initial duty cycle of the pin. Must be between 0 and 1.
- **frequency** (*float*) – The PWM frequency to use.
- **pig** (*PigpioConnection*) – pigpiod connection to use; if not specified, a new one is established

Attributes

<code>_DEFAULT_FREQUENCY</code>	<code>int([x]) -> integer</code>
<code>_DEFAULT_SOFT_FREQ</code>	<code>int([x]) -> integer</code>
<code>_HARDWARE_PWM_PINS</code>	Built-in immutable sequence.
<code>duty</code>	Returns the PWM duty cycle (pulled straight from pigpiod) mapped to the range [0-1]
<code>frequency</code>	Return the current PWM frequency active on the pin.
<code>hardware_enabled</code>	Return true if this is a hardware-enabled PWM pin; False if not.

Methods

<code>__PWMOutput__hardware_pwm(frequency, duty)</code>	Used for pins where hardware pwm is available.
<code>__PWMOutput__pwm(frequency, duty)</code>	Sets a PWM frequency and duty using either hardware or software generated PWM according to the value of <code>self.hardware_enabled</code> .

continues on next page

Table 16 – continued from previous page

<code>__PWMOutput__software_pwm(frequency, duty)</code>	Used for pins where hardware PWM is NOT available.
<code>__init__(pin[, initial_duty, frequency, pig])</code>	Inherits attributes from parent Pin, then sets PWM frequency & initial duty (use defaults if None)
<code>__duty()</code>	Returns the pigpio integer representation of the duty cycle.
<code>read()</code>	Overridden to return duty cycle instead of reading the value on the pin.
<code>write(value)</code>	Overridden to write duty cycle.

`__DEFAULT_FREQUENCY = 20000`

`__DEFAULT_SOFT_FREQ = 2000`

`__HARDWARE_PWM_PINS = (12, 13, 18, 19)`

`__init__(pin, initial_duty=0, frequency=None, pig=None)`

Inherits attributes from parent Pin, then sets PWM frequency & initial duty (use defaults if None)

Parameters

- **pin** (*int*) – The number of the pin to use. Hardware PWM pins are 12, 13, 18, and 19.
- **initial_duty** (*float*) – The initial duty cycle of the pin. Must be between 0 and 1.
- **frequency** (*float*) – The PWM frequency to use.
- **pig** (*PigpioConnection*) – pigpiod connection to use; if not specified, a new one is established

property `hardware_enabled`

Return true if this is a hardware-enabled PWM pin; False if not. The Raspberry Pi only supports hardware-generated PWM on pins 12, 13, 18, and 19, so generally `hardware_enabled` will be true if this is one of those, and false if it is not. However, `hardware_enabled` can also be dynamically set to False if for some reason pigpio is unable to start a hardware PWM (i.e. if the clock is unavailable or in use or something)

property `frequency`

Return the current PWM frequency active on the pin.

`__duty() → int`

Returns the pigpio integer representation of the duty cycle.

property `duty`

Returns the PWM duty cycle (pulled straight from pigpiod) mapped to the range [0-1]

`read() → float`

Overridden to return duty cycle instead of reading the value on the pin.

`write(value)`

Overridden to write duty cycle.

Parameters **value** (*float*) – See `PWMOutput.duty`

`__PWMOutput__hardware_pwm(frequency, duty)`

Used for pins where hardware pwm is available. -Tries to write a hardware pwm. result == 0 if it succeeds.
-Sets `hardware_enabled` flag to indicate success or failure

Parameters

- **frequency** – A new PWM frequency to use.

- **duty** (*int*) – The PWM duty cycle to set. Must be between 0 and 1.

__PWMOutput__pwm (*frequency*, *duty*)

Sets a PWM frequency and duty using either hardware or software generated PWM according to the value of *self.hardware_enabled*. If *hardware_enabled*, starts a hardware pwm with the requested duty. If not *hardware_enabled*, or if there is a problem setting a hardware generated PWM, starts a software PWM.

Parameters

- **frequency** (*float*) – A new PWM frequency to use.
- **duty** (*float*) – The PWM duty cycle to set. Must be between 0 and 1.

__PWMOutput__software_pwm (*frequency*, *duty*)

Used for pins where hardware PWM is NOT available.

Parameters

- **frequency** – A new PWM frequency to use.
- **duty** (*int*) – A pigpio integer representation of duty cycle

7.1.1.4 vent.io.devices.sensors module

Classes

<i>AnalogSensor</i> (adc, **kwargs)	Generalized class describing an analog sensor attached to the ADS1115 analog to digital converter.
<i>DLiteSensor</i> (adc, **kwargs)	D-Lite flow sensor setup.
<i>SFM3200</i> ([address, i2c_bus, pig])	I2C Inspiratory flow sensor manufactured by Sensirion AG.
<i>Sensor</i> ()	Abstract base Class describing generalized sensors.
<i>SimSensor</i> ([low, high, pig])	TODO Stub simulated sensor.

class vent.io.devices.sensors.**Sensor**

Bases: *abc.ABC*

Abstract base Class describing generalized sensors. Defines a mechanism for limited internal storage of recent observations and methods to pull that data out for external use.

Upon creation, calls *update()* to ensure that if *get* is called there will be something to return. **Attributes**

<i>__DEFAULT_STORED_OBSERVATIONS</i>	<i>int</i> ([x]) -> integer
<i>data</i>	Returns all Locally-stored observations.
<i>maxlen_data</i>	Returns the number of observations kept in the Sensor's internal ndarray.

Methods

<i>__init__</i> ()	Upon creation, calls <i>update()</i> to ensure that if <i>get</i> is called there will be something to return.
<i>_clear</i> ()	Resets the sensors internal memory.
<i>_convert</i> (raw)	Converts a raw reading from a sensor in whatever format the device communicates with into a meaningful result.

continues on next page

Table 19 – continued from previous page

<code>_raw_read()</code>	Requests a new observation from the device and returns the raw result in whatever format/units the device communicates with.
<code>_read()</code>	Calls <code>_raw_read</code> and scales the result before returning it.
<code>_verify(value)</code>	Validate reading and throw exception/alarm if sensor does not appear to be working correctly.
<code>age()</code>	Returns the time in seconds since the last sensor update, or -1 if never updated.
<code>get([average])</code>	Return the most recent sensor reading, or an average of readings since last <code>get()</code> .
<code>reset()</code>	Resets the sensors internal memory.
<code>update()</code>	Make a sensor reading, verify that it makes sense and store the result internally.

`_DEFAULT_STORED_OBSERVATIONS = 128`

`__init__()`

Upon creation, calls `update()` to ensure that if `get` is called there will be something to return.

`update()` → float

Make a sensor reading, verify that it makes sense and store the result internally. Returns True if reading was verified and False if something went wrong.

`get(average=False)` → float

Return the most recent sensor reading, or an average of readings since last `get()`. Clears internal memory so as not to have stale data.

`age()` → float

Returns the time in seconds since the last sensor update, or -1 if never updated.

`reset()`

Resets the sensors internal memory. May be overloaded by subclasses to extend functionality specific to a device.

`_clear()`

Resets the sensors internal memory.

property `data`

Returns all Locally-stored observations.

Returns An array of timestamped observations arranged oldest to newest.

Return type np.array

property `maxlen_data`

Returns the number of observations kept in the Sensor's internal ndarray. Once the ndarray has been filled, the sensor begins overwriting the oldest elements of the ndarray with new observations such that the size of the internal storage stays constant.

`_read()` → float

Calls `_raw_read` and scales the result before returning it.

abstract `_verify(value)`

Validate reading and throw exception/alarm if sensor does not appear to be working correctly.

abstract `_convert(raw)`

Converts a raw reading from a sensor in whatever format the device communicates with into a meaningful result.

```
abstract _raw_read()
```

Requests a new observation from the device and returns the raw result in whatever format/units the device communicates with.

```
_abc_impl = <_abc_data object>
```

```
class vent.io.devices.sensors.AnalogSensor(adc, **kwargs)
```

Bases: `vent.io.devices.sensors.Sensor`

Generalized class describing an analog sensor attached to the ADS1115 analog to digital converter. Inherits from the sensor base class and extends with functionality specific to analog sensors attached to the ADS1115. If instantiated without a subclass, conceptually represents a voltmeter with a normalized output.

Links analog sensor on the ADC with configuration options specified. If no options are specified, it assumes the settings currently on the ADC.

Parameters

- **adc** (`vent.io.devices.ADS1115`) – The adc object to which the AnalogSensor is attached
- ****kwargs** – *field=value* - see `vent.io.devices.ADS1115` for additional documentation. Strongly suggested to specify `MUX=adc_pin` here unless you know what you're doing.

Attributes

<code>_DEFAULT_CALIBRATION</code>	<code>dict()</code> -> new empty dictionary
<code>_DEFAULT_offset_voltage</code>	<code>int([x])</code> -> integer
<code>_DEFAULT_output_span</code>	<code>int([x])</code> -> integer

Methods

<code>__init__(adc, **kwargs)</code>	Links analog sensor on the ADC with configuration options specified.
<code>_check_and_set_attr(**kwargs)</code>	Checks to see if arguments passed to <code>__init__</code> are recognized as user configurable or calibration fields.
<code>_convert(raw)</code>	Scales raw voltage into the range 0 - 1.
<code>_fill_attr()</code>	Examines self to see if there are any fields identified as user configurable or calibration that have not been write (i.e.
<code>_raw_read()</code>	Builds kwargs from configured fields to pass along to adc, then calls <code>adc.read_conversion()</code> , which returns a raw voltage.
<code>_read()</code>	Returns a value in the range of 0 - 1 corresponding to a fraction of the full input range of the sensor.
<code>_verify(value)</code>	Checks to make sure sensor reading was indeed in [0, 1].
<code>calibrate(**kwargs)</code>	Sets the calibration of the sensor, either to the values contained in the passed tuple or by some routine; the current routine is pretty rudimentary and only calibrates offset voltage.

```
_DEFAULT_offset_voltage = 0
```

```
_DEFAULT_output_span = 5
```

```
_DEFAULT_CALIBRATION = {'conversion_factor': 1, 'offset_voltage': 0, 'output_span':
```

`__init__(adc, **kwargs)`

Links analog sensor on the ADC with configuration options specified. If no options are specified, it assumes the settings currently on the ADC.

Parameters

- **adc** (`vent.io.devices.ADS1115`) – The adc object to which the AnalogSensor is attached
- ****kwargs** – *field=value* - see `vent.io.devices.ADS1115` for additional documentation. Strongly suggested to specify *MUX=adc_pin* here unless you know what you're doing.

`calibrate(**kwargs)`

Sets the calibration of the sensor, either to the values contained in the passed tuple or by some routine; the current routine is pretty rudimentary and only calibrates offset voltage.

Parameters ****kwargs** – *calibration_field=value*, where calibration field is one of the following: 'offset_voltage', 'output_span' or 'conversion_factor'

`_read()` → `float`

Returns a value in the range of 0 - 1 corresponding to a fraction of the full input range of the sensor.

`_verify(value)` → `bool`

Checks to make sure sensor reading was indeed in [0, 1].

Parameters **value** (`float`) – Sensor reading to validate

`_convert(raw)` → `float`

Scales raw voltage into the range 0 - 1.

Parameters **raw** (`float`) – The raw sensor reading to convert.

`_raw_read()` → `float`

Builds kwargs from configured fields to pass along to adc, then calls `adc.read_conversion()`, which returns a raw voltage.

`_fill_attr()`

Examines self to see if there are any fields identified as user configurable or calibration that have not been write (i.e. were not passed to `__init__` as ****kwargs**). If a field is missing, grabs the default value either from the ADC or from `_DEFAULT_CALIBRATION` and sets it as an attribute.

`_check_and_set_attr(**kwargs)`

Checks to see if arguments passed to `__init__` are recognized as user configurable or calibration fields. If so, write the value as an attribute like: `self.KEY = VALUE`. Keeps track of how many attributes are write in this way; if at the end there unknown arguments leftover, raises a `TypeError`; otherwise, calls `_fill_attr()` to fill in fields that were not passed as arguments.

Parameters ****kwargs** – *field=value* - see `vent.io.devices.ADS1115` for additional documentation

`_abc_impl = <_abc_data object>`

class `vent.io.devices.sensors.DLiteSensor(adc, **kwargs)`

Bases: `vent.io.devices.sensors.AnalogSensor`

D-Lite flow sensor setup. This consists of the GE D-Lite sensor configured with vacuum lines running to an analog differential pressure sensor.

Links analog sensor on the ADC with configuration options specified. If no options are specified, it assumes the settings currently on the ADC.

Parameters

- **adc** (*vent.io.devices.ADS1115*) – The adc object to which the AnalogSensor is attached
- ****kwargs** – *field=value* - see *vent.io.devices.ADS1115* for additional documentation. Strongly suggested to specify *MUX=adc_pin* here unless you know what you’re doing.

Methods

<code>_convert(raw)</code>	Converts the raw differential voltage signal to a measurement of flow in liters-per-minute (LPM).
<code>calibrate(**kwargs)</code>	Do not run a calibration routine.

`_convert (raw) → float`

Converts the raw differential voltage signal to a measurement of flow in liters-per-minute (LPM).

We calibrate the D-Lite flow readings using the (pre-calibrated) Sensirion flow sensor (see SFM3200).

Parameters **raw** (*float*) – The raw sensor reading to convert.

calibrate (***kwargs*)

Do not run a calibration routine. Overrides attempt to calibrate.

`_abc_impl = <_abc_data object>`

class *vent.io.devices.sensors.SFM3200* (*address=64, i2c_bus=1, pig=None*)

Bases: *vent.io.devices.sensors.Sensor, vent.io.devices.base.I2CDevice*

I2C Inspiratory flow sensor manufactured by Sensirion AG. Range: +/- 250 SLM Datasheet:

Attributes

<code>_DEFAULT_ADDRESS</code>	<code>int([x]) -> integer</code>
<code>_FLOW_OFFSET</code>	<code>int([x]) -> integer</code>
<code>_FLOW_SCALE_FACTOR</code>	<code>int([x]) -> integer</code>

Methods

<code>__init__([address, i2c_bus, pig])</code>	param address The I2C Address of the SFM3200 (usually 0x40)
<code>_convert(raw)</code>	Overloaded to replace with device-specific protocol.
<code>_raw_read()</code>	Performs an read on the sensor, converts received bytearray, discards the last two bytes (crc values - could implement in future), and returns a signed int converted from the big endian two complement that remains.
<code>_start()</code>	Device specific:Sends the ‘start measurement’ command to the sensor.
<code>_verify(value)</code>	No further verification needed for this sensor.
<code>reset()</code>	Extended to add device specific behavior: Asks the sensor to perform a soft reset.

[https://www.sensirion.com/fileadmin/user_upload/customers/sensirion/Dokumente/ ...](https://www.sensirion.com/fileadmin/user_upload/customers/sensirion/Dokumente/.../5_Mass_Flow_Meters/Datasheets/Sensirion_Mass_Flow_Meters_SFM3200_Datasheet.pdf)
[5_Mass_Flow_Meters/Datasheets/Sensirion_Mass_Flow_Meters_SFM3200_Datasheet.pdf](https://www.sensirion.com/fileadmin/user_upload/customers/sensirion/Dokumente/.../5_Mass_Flow_Meters/Datasheets/Sensirion_Mass_Flow_Meters_SFM3200_Datasheet.pdf)

Parameters

- **address** (*int*) – The I2C Address of the SFM3200 (usually 0x40)
- **i2c_bus** (*int*) – The I2C Bus to use (usually 1 on the Raspberry Pi)
- **pig** (*PigpioConnection*) – pigpiod connection to use; if not specified, a new one is established

```
_DEFAULT_ADDRESS = 64
```

```
_FLOW_OFFSET = 32768
```

```
_FLOW_SCALE_FACTOR = 120
```

```
__init__ (address=64, i2c_bus=1, pig=None)
```

Parameters

- **address** (*int*) – The I2C Address of the SFM3200 (usually 0x40)
- **i2c_bus** (*int*) – The I2C Bus to use (usually 1 on the Raspberry Pi)
- **pig** (*PigpioConnection*) – pigpiod connection to use; if not specified, a new one is established

```
reset ()
```

Extended to add device specific behavior: Asks the sensor to perform a soft reset. 80 ms soft reset time.

```
_start ()
```

Device specific: Sends the ‘start measurement’ command to the sensor. Start-up time once command has been recieved is ‘less than 100ms’

```
_verify (value) → bool
```

No further verification needed for this sensor. Onboard chip handles all that. Could throw in a CRC8 checker instead of discarding them in `_convert()`.

Parameters **value** (*float*) – The sensor reading to verify.

```
_convert (raw) → float
```

Overloaded to replace with device-specific protocol. Convert raw int to a flow reading having type float with units slm. Implementation differs from parent for clarity and consistency with source material.

Source:

https://www.sensirion.com/fileadmin/user_upload/customers/sensirion/Dokumente/...

5_Mass_Flow_Meters/Application_Notes/Sensirion_Mass_Flo_Meters_SFM3xxx_I2C_Functional_Description.

Parameters **raw** (*int*) – The raw value read from the SFM3200

```
_raw_read () → int
```

Performs an read on the sensor, converts received bytearray, discards the last two bytes (crc values - could implement in future), and returns a signed int converted from the big endian two complement that remains.

```
_abc_impl = <_abc_data object>
```

```
class vent.io.devices.sensors.SimSensor (low=0, high=100, pig=None)
```

Bases: `vent.io.devices.sensors.Sensor`

TODO Stub simulated sensor.

Parameters

- **low** – Lower-bound of possible sensor values

- **high** – Upper-bound of possible sensor values
- **pig** (*PigpioConnection*) – Ignored.

Methods

<code>__init__([low, high, pig])</code>	TODO Stub simulated sensor.
<code>_convert(raw)</code>	Does nothing for a simulated sensor.
<code>_raw_read()</code>	Initializes randomly, otherwise does a random walk-ish thing.
<code>_verify(value)</code>	Usually verifies sensor readings but occasionally misbehaves.

`__init__ (low=0, high=100, pig=None)`
 TODO Stub simulated sensor.

Parameters

- **low** – Lower-bound of possible sensor values
- **high** – Upper-bound of possible sensor values
- **pig** (*PigpioConnection*) – Ignored.

`_verify (value) → bool`

Usually verifies sensor readings but occasionally misbehaves.

Parameters **value** (*float*) – The sensor reading to verify

`_convert (raw) → float`

Does nothing for a simulated sensor. Returns what it is passed.

Parameters **raw** (*float*) – The raw value to convert

`_abc_impl = <_abc_data object>`

`_raw_read () → float`

Initializes randomly, otherwise does a random walk-ish thing.

7.1.1.5 vent.io.devices.valves module

Classes

<code>OnOffValve(pin[, form, pig])</code>	An extension of <code>vent.io.iobase.Pin</code> which uses valve terminology for its methods.
<code>PWMControlValve(pin[, form, frequency, ...])</code>	An extension of <code>PWMOutput</code> which incorporates linear compensation of the valve's response.
<code>SimControlValve([pin, form, frequency, ...])</code>	stub: a simulated linear control valve
<code>SimOnOffValve([pin, form, pig])</code>	stub: a simulated on/off valve
<code>SolenoidBase([form])</code>	An abstract baseclass that defines methods using valve terminology.

class `vent.io.devices.valves.SolenoidBase (form='Normally Closed')`

Bases: `abc.ABC`

An abstract baseclass that defines methods using valve terminology. Also allows configuring both normally open and normally closed valves (called the “form” of the valve).

Parameters **form** (*str*) – The form of the solenoid; can be either *Normally Open* or *Normally*

Closed

Attributes

<code>__FORMS</code>	dict() -> new empty dictionary
<code>form</code>	Returns the human-readable form of the valve.
<code>is_open</code>	Returns True if valve is open, False if it is closed

Methods

<code>__init__([form])</code>	param form The form of the solenoid; can be either <i>Normally Open</i> or <i>Normally Closed</i>
<code>close()</code>	De-energizes valve if Normally Closed.
<code>open()</code>	Energizes valve if Normally Closed.

```
__FORMS = {'Normally Closed': 0, 'Normally Open': 1}
```

```
__init__(form='Normally Closed')
```

Parameters form (*str*) – The form of the solenoid; can be either *Normally Open* or *Normally Closed*

property form

Returns the human-readable form of the valve.

abstract open()

Energizes valve if Normally Closed. De-energizes if Normally Open.

abstract close()

De-energizes valve if Normally Closed. Energizes if Normally Open.

abstract property is_open

Returns True if valve is open, False if it is closed

```
__abc_impl = <_abc_data object>
```

```
class vent.io.devices.valves.OnOffValve(pin, form='Normally Closed', pig=None)
```

Bases: `vent.io.devices.valves.SolenoidBase`, `vent.io.devices.pins.Pin`

An extension of `vent.io.iobase.Pin` which uses valve terminology for its methods. Also allows configuring both normally `_open` and normally closed valves (called the “form” of the valve).

Parameters

- **pin** (*int*) – The number of the pin to use
- **form** (*str*) – The form of the solenoid; can be either *Normally Open* or *Normally Closed*
- **pig** (*PigpioConnection*) – pigpiod connection to use; if not specified, a new one is established

Attributes

<code>__FORMS</code>	dict() -> new empty dictionary
----------------------	--------------------------------

continues on next page

Table 29 – continued from previous page

<code>is_open</code>	Implements parent's abstractmethod; returns True if valve is open, False if it is closed
----------------------	--

Methods

<code>__init__(pin[, form, pig])</code>	param pin The number of the pin to use
<code>close()</code>	De-energizes valve if Normally Closed.
<code>open()</code>	Energizes valve if Normally Closed.

```
_FORMS = {'Normally Closed': 0, 'Normally Open': 1}
```

```
__init__ (pin, form='Normally Closed', pig=None)
```

Parameters

- **pin** (*int*) – The number of the pin to use
- **form** (*str*) – The form of the solenoid; can be either *Normally Open* or *Normally Closed*
- **pig** (*PigpioConnection*) – pigpiod connection to use; if not specified, a new one is established

open()
Energizes valve if Normally Closed. De-energizes if Normally Open.

close()
De-energizes valve if Normally Closed. Energizes if Normally Open.

property is_open
Implements parent's abstractmethod; returns True if valve is open, False if it is closed

```
_abc_impl = <_abc_data object>
```

```
class vent.io.devices.valves.PWMControlValve (pin, form='Normally Closed', frequency=None, response=None, pig=None)
```

Bases: `vent.io.devices.valves.SolenoidBase`, `vent.io.devices.pins.PWMOutput`

An extension of PWMOutput which incorporates linear compensation of the valve's response.

Parameters

- **pin** (*int*) – The number of the pin to use
- **form** (*str*) – The form of the solenoid; can be either *Normally Open* or *Normally Closed*
- **frequency** (*float*) – The PWM frequency to use.
- **response** (*str*) – “/path/to/response/curve/file”
- **pig** (*PigpioConnection*) – pigpiod connection to use; if not specified, a new one is established

Methods

<code>__init__(pin[, form, frequency, response, pig])</code>	param pin The number of the pin to use
<code>_load_valve_response(response_path)</code>	Loads and applies a response curve of the form $f(\text{setpoint}) = \text{duty}$.
<code>close()</code>	Implements parent's abstractmethod; fully closes the valve
<code>inverse_response(duty_cycle[, rising])</code>	Inverse of response.
<code>open()</code>	Implements parent's abstractmethod; fully opens the valve
<code>response(setpoint[, rising])</code>	Setpoint takes a value in the range (0,100) so as not to confuse with duty cycle, which takes a value in the range (0,1).

Attributes

<code>is_open</code>	Implements parent's abstractmethod; returns True if valve is open, False if it is closed
<code>setpoint</code>	The linearized setpoint corresponding to the current duty cycle according to the valve's response curve

`__init__(pin, form='Normally Closed', frequency=None, response=None, pig=None)`

Parameters

- **pin** (*int*) – The number of the pin to use
- **form** (*str*) – The form of the solenoid; can be either *Normally Open* or *Normally Closed*
- **frequency** (*float*) – The PWM frequency to use.
- **response** (*str*) – “/path/to/response/curve/file”
- **pig** (*PigpioConnection*) – pigpiod connection to use; if not specified, a new one is established

property is_open

Implements parent's abstractmethod; returns True if valve is open, False if it is closed

open()

Implements parent's abstractmethod; fully opens the valve

close()

Implements parent's abstractmethod; fully closes the valve

property setpoint

The linearized setpoint corresponding to the current duty cycle according to the valve's response curve

Returns A number between 0 and 1 representing the current flow as a proportion of maximum

Return type *float*

response(setpoint, rising=True)

Setpoint takes a value in the range (0,100) so as not to confuse with duty cycle, which takes a value in the range (0,1). Response curves are specific to individual valves and are to be implemented by subclasses. Different curves are calibrated to 'rising = True' (valves opening) or 'rising = False' (valves closing), as different characteristic flow behavior can be observed.

Parameters

- **setpoint** (*float*) – A number between 0 and 1 representing how much to open the valve
- **rising** (*bool*) – Whether or not the requested setpoint is higher than the last (rising = True), or the opposite (Rising = False)

Returns The PWM duty cycle corresponding to the requested setpoint

Return type *float*

inverse_response (*duty_cycle*, *rising=True*)

Inverse of response. Given a duty cycle in the range (0,1), returns the corresponding linear setpoint in the range (0,100).

Parameters

- **duty_cycle** – The PWM duty cycle
- **rising** (*bool*) – Whether or not the requested setpoint is higher than the last (rising = True), or the opposite (Rising = False)

Returns The setpoint of the valve corresponding to *duty_cycle*

Return type *float*

_load_valve_response (*response_path*)

Loads and applies a response curve of the form $f(\text{setpoint}) = \text{duty}$. A response curve maps the underlying PWM duty cycle *duty* onto the normalized variable *setpoint* representing the flow through the valve as a percentage of its maximum.

Flow through a proportional valve may be nonlinear with respect to [PWM] duty cycle, if the valve itself does not include its own electronics to linearize response wrt/ input. Absent on-board compensation of response, a proportional solenoid with likely not respond [flow] at all below some minimum threshold duty cycle. Above this threshold, the proportional valve begins to open and its response resembles a sigmoid: just past the threshold there is a region where flow increases exponentially wrt/ duty cycle, this is followed by a region of pseudo-linear response that begins to taper off, eventually approaching the valve's maximum flow asymptotically as the duty cycle approaches 100% and the valve opens fully.

Parameters **response_path** – 'path/to/binary/response/file' - if response_path is None, defaults to *setpoint = duty*

_abc_impl = <_abc_data object>

class `vent.io.devices.valves.SimOnOffValve` (*pin=None*, *form='Normally Closed'*, *pig=None*)

Bases: `vent.io.devices.valves.SolenoidBase`

stub: a simulated on/off valve

Args: form (str): The form of the solenoid; can be either *Normally Open* or *Normally Closed* **Attributes**

<i>is_open</i>	Returns True if valve is open, False if it is closed
----------------	--

Methods

<i>close()</i>	De-energizes valve if Normally Closed.
<i>open()</i>	Energizes valve if Normally Closed.

open ()

Energizes valve if Normally Closed. De-energizes if Normally Open.

close()

De-energizes valve if Normally Closed. Energizes if Normally Open.

property is_open

Returns True if valve is open, False if it is closed

_abc_impl = <_abc_data object>

```
class vent.io.devices.valves.SimControlValve (pin=None, form='Normally Closed',
                                             frequency=None, response=None,
                                             pig=None)
```

Bases: `vent.io.devices.valves.SolenoidBase`

stub: a simulated linear control valve

Parameters

- **pin** (*int*) – (unused for sim)
- **form** (*str*) – The form of the solenoid; can be either *Normally Open* or *Normally Closed*
- **frequency** (*float*) – (unused for sim)
- **response** (*str*) – (unused for sim) # TODO implement this (requires refactor)
- **pig** (*PigpioConnection*) – (unused for sim)

Methods

`__init__`([pin, form, frequency, response, pig])

param pin (unused for sim)

<code>close()</code>	Implements parent's abstractmethod; fully closes the valve
----------------------	--

<code>open()</code>	Implements parent's abstractmethod; fully opens the valve
---------------------	---

Attributes

<code>is_open</code>	Implements parent's abstractmethod; returns True if valve is open, False if it is closed
----------------------	--

<code>setpoint</code>	The requested linearized set-point of the valve.
-----------------------	--

`__init__` (pin=None, form='Normally Closed', frequency=None, response=None, pig=None)

Parameters

- **pin** (*int*) – (unused for sim)
- **form** (*str*) – The form of the solenoid; can be either *Normally Open* or *Normally Closed*
- **frequency** (*float*) – (unused for sim)
- **response** (*str*) – (unused for sim) # TODO implement this (requires refactor)
- **pig** (*PigpioConnection*) – (unused for sim)

property is_open

Implements parent's abstractmethod; returns True if valve is open, False if it is closed FIXME: Needs refactor; duplicate property to PWMControlValve.is_open

_abc_impl = <_abc_data object>

open()

Implements parent's abstractmethod; fully opens the valve FIXME: Needs refactor; duplicate method to PWMControlValve.open()

close()

Implements parent's abstractmethod; fully closes the valve FIXME: Needs refactor; duplicate method to PWMControlValve.close()

property setpoint

The requested linearized set-point of the valve.

Returns A number between 0 and 1 representing the current flow as a proportion of maximum

Return type float

7.1.1.6 Module contents

A module for ventilator hardware device drivers

7.2 Submodules

7.3 vent.io.hal module

Module for interacting with physical and/or simulated devices installed on the ventilator.

Classes

<code>Hal([config_file])</code>	Hardware Abstraction Layer for ventilator hardware.
---------------------------------	---

class `vent.io.hal.Hal` (*config_file*='vent/io/config/devices.ini')

Bases: `object`

Hardware Abstraction Layer for ventilator hardware. Defines a common API for interacting with the sensors & actuators on the ventilator. The types of devices installed on the ventilator (real or simulated) are specified in a configuration file.

Initializes HAL from config_file. For each section in *config_file*, imports the class <type> from module <module>, and sets attribute `self.<section>` = <type>(***opts*), where *opts* is a dict containing all of the options in <section> that are not <type> or <section>. For example, upon encountering the following entry in *config_file.ini*:

```
[adc] type = ADS1115 module = devices i2c_address = 0x48 i2c_bus = 1
```

The Hal will:

- 1) **Import vent.io.devices.ADS1115 (or ADS1015) as a local variable:** `class_ = getattr(import_module('devices', 'vent.io'), 'ADS1115')`
- 2) **Instantiate an ADS1115 object with the arguments defined in config_file and set it as an attribute:** `self._adc = class_(pig=self._pig,address=0x48,i2c_bus=1)`

Note: `RawConfigParser.optionxform()` is overloaded here s.t. options are case sensitive (they are by default case insensitive). This is necessary due to the kwarg `MUX` which is so named for consistency with the config registry documentation in the ADS1115 datasheet. For example, A P4vMini pressure_sensor on pin A0 (`MUX=0`) of the ADC is passed arguments like:

```
analog_sensor = AnalogSensor( pig=self._pig, adc=self._adc, MUX=0, offset_voltage=0.25, out-  
put_span = 4.0, conversion_factor=2.54*20  
)
```

Note: `ast.literal_eval(opt)` interprets integers, `0xFF`, `(a, b)` etc. correctly. It does not interpret strings correctly, nor does it know `'adc' -> self._adc`; therefore, these special cases are explicitly handled.

Methods

<code>__init__([config_file])</code>	Initializes HAL from <code>config_file</code> .
--------------------------------------	---

Attributes

<code>aux_pressure</code>	Returns the pressure from the auxiliary pressure sensor, if so equipped.
<code>flow_ex</code>	The measured flow rate expiratory side.
<code>flow_in</code>	The measured flow rate inspiratory side.
<code>pressure</code>	Returns the pressure from the primary pressure sensor.
<code>setpoint_ex</code>	The currently requested flow on the expiratory side as a proportion of the maximum.
<code>setpoint_in</code>	The currently requested flow for the inspiratory proportional control valve as a proportion of maximum.

Parameters `config_file` (*str*) – Path to the configuration file containing the definitions of specific components on the ventilator machine. (e.g., `config_file = "vent/io/config/devices.ini"`)

```
__init__(config_file='vent/io/config/devices.ini')
```

Initializes HAL from `config_file`. For each section in `config_file`, imports the class `<type>` from module `<module>`, and sets attribute `self.<section> = <type>(**opts)`, where `opts` is a dict containing all of the options in `<section>` that are not `<type>` or `<section>`. For example, upon encountering the following entry in `config_file.ini`:

```
[adc] type = ADS1115 module = devices i2c_address = 0x48 i2c_bus = 1
```

The Hal will:

- 1) **Import `vent.io.devices.ADS1115` (or `ADS1015`) as a local variable:** `class_ = getattr(import_module('devices', 'vent.io'), 'ADS1115')`
- 2) **Instantiate an `ADS1115` object with the arguments defined in `config_file` and set it as an attribute:** `self._adc = class_(pig=self._pig,address=0x48,i2c_bus=1)`

Note: `RawConfigParser.optionxform()` is overloaded here s.t. options are case sensitive (they are by default case insensitive). This is necessary due to the kwarg `MUX` which is so named for consistency with the config registry documentation in the ADS1115 datasheet. For example, A P4vMini pressure_sensor on pin A0 (`MUX=0`) of the ADC is passed arguments like:

```
analog_sensor = AnalogSensor( pig=self._pig, adc=self._adc, MUX=0, offset_voltage=0.25, out-  
    put_span = 4.0, conversion_factor=2.54*20  
    )
```

Note: `ast.literal_eval(opt)` interprets integers, 0xFF, (a, b) etc. correctly. It does not interpret strings correctly, nor does it know 'adc' -> `self._adc`; therefore, these special cases are explicitly handled.

Parameters `config_file` (*str*) – Path to the configuration file containing the definitions of specific components on the ventilator machine. (e.g., `config_file = "vent/io/config/devices.ini"`)

property pressure

Returns the pressure from the primary pressure sensor.

property aux_pressure

Returns the pressure from the auxiliary pressure sensor, if so equipped. If a secondary pressure sensor is not defined, raises a `RuntimeWarning`.

property flow_in

The measured flow rate inspiratory side.

property flow_ex

The measured flow rate expiratory side.

property setpoint_in

The currently requested flow for the inspiratory proportional control valve as a proportion of maximum.

property setpoint_ex

The currently requested flow on the expiratory side as a proportion of the maximum.

7.4 Module contents

8.1 Main Alarm Module

Classes

AlarmSeverity	An enumeration.
AlarmType	An enumeration.

class vent.alarm.**AlarmType**
An enumeration. **Attributes**

HIGH_O2	int([x]) -> integer
HIGH_PEEP	int([x]) -> integer
HIGH_PRESSURE	int([x]) -> integer
HIGH_VTE	int([x]) -> integer
LEAK	int([x]) -> integer
LOW_O2	int([x]) -> integer
LOW_PEEP	int([x]) -> integer
LOW_PRESSURE	int([x]) -> integer
LOW_VTE	int([x]) -> integer
OBSTRUCTION	int([x]) -> integer

```
LOW_PRESSURE = 1
HIGH_PRESSURE = 2
LOW_VTE = 3
HIGH_VTE = 4
LOW_PEEP = 5
HIGH_PEEP = 6
LOW_O2 = 7
HIGH_O2 = 8
OBSTRUCTION = 9
LEAK = 10
```

class vent.alarm.**AlarmSeverity**
An enumeration. **Attributes**

HIGH	int([x]) -> integer
LOW	int([x]) -> integer
MEDIUM	int([x]) -> integer
OFF	int([x]) -> integer

HIGH = 3

MEDIUM = 2

LOW = 1

OFF = 0

8.2 Alarm Manager

Classes

Alarm_Manager()	vent.alarm.alarm_manager.active_alarms
-----------------	--

class vent.alarm.alarm_manager.**Alarm_Manager**

Attributes

active_alarms	dict() -> new empty dictionary
callbacks	Built-in mutable sequence.
cleared_alarms	Built-in mutable sequence.
dependencies	dict() -> new empty dictionary
logged_alarms	Built-in mutable sequence.
pending_clears	Built-in mutable sequence.
rules	dict() -> new empty dictionary
snoozed_alarms	dict() -> new empty dictionary

Methods

add_callback(callback)	
check_rule(rule, sensor_values)	
clear_all_alarms()	
deactivate_alarm(alarm)	Mark an alarm's internal active flags and remove from active_alarms
dismiss_alarm(alarm_type, duration)	GUI or other object requests an alarm to be dismissed & deactivated
emit_alarm(alarm_type, severity)	Emit alarm (by calling all callbacks with it).
get_alarm_severity(alarm_type)	
load_rule(alarm_rule)	
load_rules()	
register_alarm(alarm)	Add alarm to registry.
register_dependency(condition, dependency)	Add dependency in a Condition object to be updated when values are changed

continues on next page

Table 6 – continued from previous page

<code>reset()</code>	reset all conditions, callbacks, and other stateful attributes and clear alarms
<code>update(sensor_values)</code>	
<code>update_dependencies(control_setting)</code>	Update Condition objects that update their value according to some control parameter

active_alarms

{AlarmType: Alarm}

Type dict

pending_clears

[AlarmType] list of alarms that have been requested to be cleared

Type list

callbackslist of callables that accept *Alarm* s when they are raised/alterd.

Type list

cleared_alarms

of AlarmType s, alarms that have been cleared but have not dropped back into the ‘off’ range to enable re-raising

Type list

snoozed_alarms

of AlarmType s : times, alarms that should not be raised because they have been silenced for a period of time

Type dict

If an Alarm_Manager already exists, when initing just return that one

_instance = None**active_alarms**: Dict[vent.alarm.AlarmType, vent.alarm.alarm.Alarm] = {}**logged_alarms**: List[vent.alarm.alarm.Alarm] = []**dependencies** = {}**pending_clears** = []**cleared_alarms** = []**snoozed_alarms** = {}**callbacks** = []**rules** = {}**load_rules**()**load_rule**(alarm_rule: vent.alarm.rule.Alarm_Rule)**update**(sensor_values: vent.common.message.SensorValues)**check_rule**(rule: vent.alarm.rule.Alarm_Rule, sensor_values: vent.common.message.SensorValues)**emit_alarm**(alarm_type: vent.alarm.AlarmType, severity: vent.alarm.AlarmSeverity)

Emit alarm (by calling all callbacks with it).

Note: This method emits *and* clears alarms – a cleared alarm is emitted with `AlarmSeverity.OFF`

Parameters

- **alarm_type** (`AlarmType`) –
- **severity** (`AlarmSeverity`) –

deactivate_alarm (*alarm*: (<enum 'AlarmType'>, <class 'vent.alarm.alarm.Alarm'>))

Mark an alarm's internal active flags and remove from `active_alarms`

Note: This does *not* alert listeners that an alarm has been cleared, for that emit an alarm with `AlarmSeverity.OFF`

Parameters alarm –

Returns:

dismiss_alarm (*alarm_type*: `vent.alarm.AlarmType`, *duration*: `float = None`)

GUI or other object requests an alarm to be dismissed & deactivated

GUI will wait until it receives an *emit_alarm* of severity == `OFF` to remove alarm widgets. If the alarm is not latched

If the alarm is latched, `alarm_manager` will not decrement alarm severity or emit *OFF* until a) the condition returns to *OFF*, and b) the user dismisses the alarm

Parameters

- **alarm_type** (`AlarmType`) – Alarm to dismiss
- **duration** (`float`) – seconds - amount of time to wait before alarm can be re-raised
If a duration is provided, the alarm will not be able to be re-raised

get_alarm_severity (*alarm_type*: `vent.alarm.AlarmType`)

register_alarm (*alarm*: `vent.alarm.alarm.Alarm`)

Add alarm to registry.

Parameters alarm (`Alarm`) –

register_dependency (*condition*: `vent.alarm.condition.Condition`, *dependency*: `dict`)

Add dependency in a `Condition` object to be updated when values are changed

Parameters

- **condition** –
- **dependency** (`dict`) – either a (`ValueName`, `attribute_name`) or optionally also + transformation callable

update_dependencies (*control_setting*: `vent.common.message.ControlSetting`)

Update `Condition` objects that update their value according to some control parameter

Parameters control_setting (`ControlSetting`) –

Returns:

add_callback (*callback*: `Callable`)

```
clear_all_alarms()
```

```
reset()
```

reset all conditions, callbacks, and other stateful attributes and clear alarms

8.3 Alarm

Classes

<code>Alarm(alarm_type, severity, start_time[, ...])</code>	Class used by the program to control and coordinate alarms.
---	---

```
class vent.alarm.alarm.Alarm(alarm_type: vent.alarm.AlarmType, severity:
                             vent.alarm.AlarmSeverity, start_time: float = None, latch=True,
                             persistent=True, value=None, message=None, managed=False)
```

Class used by the program to control and coordinate alarms.

Parameterized by a `Alarm_Rule` and managed by `Alarm_Manager`

Parameters

- **alarm_type** –
- **severity** –
- **start_time** –
- **value** (*int*, *float*) – optional - numerical value that generated the alarm
- **message** (*str*) – optional - override default text generated by `AlarmManager`
- **managed** (*bool*) – if created by the `alarm_manager`, don't register

Methods

```
__init__(alarm_type, severity, start_time[, ...])
```

param alarm_type

```
deactivate()
```

Attributes

<code>alarm_type</code>	
<code>id_counter</code>	used to generate unique IDs for each alarm
<code>manager</code>	have ta do it this janky way to avoid circular imports
<code>severity</code>	

```
id_counter = count(0)
```

used to generate unique IDs for each alarm

Type `itertools.count`

```
__init__(alarm_type: vent.alarm.AlarmType, severity: vent.alarm.AlarmSeverity, start_time: float =
         None, latch=True, persistent=True, value=None, message=None, managed=False)
```

Parameters

- **alarm_type** –

- **severity** –
- **start_time** –
- **value** (*int*, *float*) – optional - numerical value that generated the alarm
- **message** (*str*) – optional - override default text generated by AlarmManager
- **managed** (*bool*) – if created by the alarm_manager, don't register

property manager

have to do it this janky way to avoid circular imports

property severity

property alarm_type

deactivate()

8.4 Condition

Classes

AlarmSeverityCondition(alarm_type, severity, ...)	param alarm_type
Condition(depends, *args, **kwargs)	Base class for specifying alarm test conditions
CycleAlarmSeverityCondition(n_cycles, *args, ...)	alarm goes out of range for a specific number of breath cycles
CycleValueCondition(n_cycles, *args, **kwargs)	value goes out of range for a specific number of breath cycles
TimeValueCondition(time, *args, **kwargs)	value goes out of range for specific amount of time
ValueCondition(value_name, limit, mode, ...)	value is greater or lesser than some max/min

Functions

get_alarm_manager()

`vent.alarm.condition.get_alarm_manager()`

class `vent.alarm.condition.Condition` (*depends: dict = None*, *args, **kwargs)

Base class for specifying alarm test conditions

Need to be able to condition alarms based on * value ranges * value ranges & durations * levels of other alarms

Methods

<code>__init__(depends, *args, **kwargs)</code>	param depends
<code>check(sensor_values)</code>	
<code>reset()</code>	If a condition is stateful, need to provide some method of resetting the state

Attributes

manager

manager

alarm manager, used to get status of alarms

Type `vent.alarm.alarm_manager.Alarm_Manager`**_child**

if another condition is added to this one, store a reference to it

Type `Condition`**Parameters**

- **depends** (*list*, *dict*) – a list of, or a single dict:

```
{'value_name': ValueName,
 'value_attr': attr in ControlMessage,
 'condition_attr',
 optional: transformation: callable)
that declare what values are needed to update
```

- ***args** –
- ****kwargs** –

__init__ (*depends*: *dict* = *None*, **args*, ***kwargs*)**Parameters**

- **depends** (*list*, *dict*) – a list of, or a single dict:

```
{'value_name': ValueName,
 'value_attr': attr in ControlMessage,
 'condition_attr',
 optional: transformation: callable)
that declare what values are needed to update
```

- ***args** –
- ****kwargs** –

property manager**check** (*sensor_values*)**reset** ()

If a condition is stateful, need to provide some method of resetting the state

```
class vent.alarm.condition.ValueCondition (value_name: vent.common.values.ValueName,
                                           limit: (<class 'int'>, <class 'float'>), mode: str,
                                           *args, **kwargs)
```

value is greater or lesser than some max/min

Parameters

- **value_name** (*ValueName*) – Which value to check
- **limit** (*int*, *float*) – value to check against
- **mode** (*'min'*, *'max'*) – whether the limit is a minimum or maximum
- ***args** –

- ****kwargs** –

Methods

<code>__init__(value_name, limit, mode, *args, ...)</code>	param value_name Which value to check
--	--

<code>check(sensor_values)</code>

<code>reset()</code>	not stateful, do nothing.
----------------------	---------------------------

Attributes

<code>mode</code>

`__init__` (*value_name*: *vent.common.values.ValueName*, *limit*: (<class 'int'>, <class 'float'>), *mode*: *str*, **args*, ***kwargs*)

Parameters

- **value_name** (*ValueName*) – Which value to check
- **limit** (*int*, *float*) – value to check against
- **mode** ('min', 'max') – whether the limit is a minimum or maximum
- ***args** –
- ****kwargs** –

property mode

check (*sensor_values*)

reset ()
not stateful, do nothing.

class *vent.alarm.condition.CycleValueCondition* (*n_cycles*, **args*, ***kwargs*)
value goes out of range for a specific number of breath cycles **Methods**

<code>check(sensor_values)</code>	
<code>reset()</code>	not stateful, do nothing.

Attributes

<code>n_cycles</code>

_start_cycle
The breath cycle where the

Type *int*

_mid_check
whether a value has left the acceptable range and we are counting consecutive breath cycles

Type *bool*

Args: *value_name* (*ValueName*): Which value to check *limit* (*int*, *float*): value to check against *mode* ('min', 'max'): whether the limit is a minimum or maximum **args*: ***kwargs*:

property `n_cycles`

check (*sensor_values*)

reset ()

not stateful, do nothing.

class `vent.alarm.condition.TimeValueCondition` (*time*, **args*, ***kwargs*)

value goes out of range for specific amount of time

Parameters

- **time** (*float*) – number of seconds value must be out of range
- ***args** –
- ****kwargs** –

Methods

`__init__` (*time*, **args*, ***kwargs*)

param time number of seconds value must be out of range

`check` (*sensor_values*)

`reset` () not stateful, do nothing.

`__init__` (*time*, **args*, ***kwargs*)

Parameters

- **time** (*float*) – number of seconds value must be out of range
- ***args** –
- ****kwargs** –

check (*sensor_values*)

reset ()

not stateful, do nothing.

class `vent.alarm.condition.AlarmSeverityCondition` (*alarm_type*: `vent.alarm.AlarmType`,
severity: `vent.alarm.AlarmSeverity`,
mode: *str* = 'min', **args*, ***kwargs*)

Parameters

- **alarm_type** –
- **severity** –
- **mode** (*str*) – one of 'min', 'equals', or 'max'. 'min' returns true if the alarm is at least this value (note the difference from ValueCondition which returns true if the alarm is less than..) and vice versa for 'max'.

Note: 'min' and 'max' use >= and <= rather than > and <

- ***args** –
- ****kwargs** –

Methods

`__init__(alarm_type, severity, mode, *args, ...)`**param alarm_type**

`check(sensor_values)``reset()`

If a condition is stateful, need to provide some method of resetting the state

Attributes

`mode`

```
__init__(alarm_type: vent.alarm.AlarmType, severity: vent.alarm.AlarmSeverity, mode: str = 'min',
          *args, **kwargs)
```

Parameters

- **alarm_type** –
- **severity** –
- **mode** (*str*) – one of ‘min’, ‘equals’, or ‘max’. ‘min’ returns true if the alarm is at least this value (note the difference from ValueCondition which returns true if the alarm is less than..) and vice versa for ‘max’.

Note: ‘min’ and ‘max’ use `>=` and `<=` rather than `>` and `<`

- ***args** –
- ****kwargs** –

property mode

`check(sensor_values)``reset()`

If a condition is stateful, need to provide some method of resetting the state

```
class vent.alarm.condition.CycleAlarmSeverityCondition(n_cycles, *args, **kwargs)
    alarm goes out of range for a specific number of breath cycles
```

Todo: note that this is exactly the same as CycleValueCondition. Need to do the multiple inheritance thing

Methods

`check(sensor_values)``reset()`

If a condition is stateful, need to provide some method of resetting the state

Attributes

`n_cycles`

_start_cycle

The breath cycle where the

Type `int`

_mid_check

whether a value has left the acceptable range and we are counting consecutive breath cycles

Type `bool`

Args: alarm_type: severity: mode (str): one of 'min', 'equals', or 'max'.

'min' returns true if the alarm is at least this value (note the difference from ValueCondition which returns true if the alarm is less than..) and vice versa for 'max'.

Note: 'min' and 'max' use `>=` and `<=` rather than `>` and `<`

***args: **kwargs:**

property `n_cycles`

check (*sensor_values*)

reset ()

If a condition is stateful, need to provide some method of resetting the state

8.5 Alarm Rule

Class to declare alarm rules

Classes

`Alarm_Rule(name, conditions[, latch, ...])`

- name of rule
-

class `vent.alarm.rule.Alarm_Rule` (*name: vent.alarm.AlarmType, conditions, latch=True, persistent=True, technical=False*)

- name of rule
- conditions: ((alarm_type, (condition_1, condition_2)), ...)
- persistent (bool): if True, alarm will not be visually dismissed until alarm conditions are no longer true
- latch (bool): if True, alarm severity cannot be decremented until user manually dismisses
- silencing/overriding rules

Methods

`check(sensor_values)`

Check all of our conditions .

`reset()`

Attributes

severity

Last Alarm Severity from .check ()

check (*sensor_values*)

Check all of our conditions .

Parameters **sensor_values** –

Returns:

property **severity**

Last Alarm Severity from .check () :returns: AlarmSeverity

reset ()

REQUIREMENTS

DATASHEETS & MANUALS

10.1 Manuals

- [Hamilton T1 Quick Guide](#)

10.2 Other Reference Material

- [Hamilton UI Simulator](#)

CHAPTER
ELEVEN

SPECS

CHANGELOG

12.1 Version 0.0

12.1.1 v0.0.2 (April xxth, 2020)

- Refactored gui into a module, splitting widgets, styles, and defaults.

12.1.2 v0.0.1 (April 12th, 2020)

- Added changelog
- Moved requirements for building docs to *requirements_docs.txt* so regular program reqs are a bit lighter.
- added autosummaries
- added additional resources & documentation files, with examples for adding external files like pdfs

12.1.3 v0.0.0 (April 12th, 2020)

Example of a changelog entry!!!

- We fixed this
- and this
- and this

Warning: but we didn't do this thing

Todo: and we still have to do this other thing.

BUILDING THE DOCS

A very brief summary...

- Docs are configured to be built from `_docs` into `docs`.
- The main page is `index.rst` which links to the existing modules
- To add a new page, you can create a new `.rst` file if you are writing with [Restructuredtext](#), or a `.md` file if you are writing with markdown.

13.1 Local Build

- `pip install -r requirements.txt`
- `cd _docs`
- `make html`

Documentation will be generated into `docs`

Advertisement :)

- [pica](#) - high quality and fast image resize in browser.
- [babelfish](#) - developer friendly i18n with plurals support and easy syntax.

You will like those projects!

H1 HEADING 8-)

14.1 h2 Heading


14.1.1 h3 Heading

14.1.1.1 h4 Heading

h5 Heading

h6 Heading

14.2 Horizontal Rules



14.3 Emphasis

This is bold text

This is bold text

This is italic text

This is italic text

14.4 Blockquotes

Blockquotes can also be nested. . .

... by using additional greater-than signs right next to each other. . .

... or with spaces between arrows.

14.5 Lists

Unordered

- Create a list by starting a line with +, −, or *
- Sub-lists are made by indenting 2 spaces:
 - Marker character change forces new list start:
 - * Ac tristique libero volutpat at
 - * Facilisis in pretium nisl aliquet
 - * Nulla volutpat aliquam velit
- Very easy!

Ordered

1. Lorem ipsum dolor sit amet
2. Consectetur adipiscing elit
3. Integer molestie lorem at massa
4. You can use sequential numbers...
5. ...or keep all the numbers as 1.

14.6 Code

Inline code

Indented code

```
// Some comments  
line 1 of code  
line 2 of code  
line 3 of code
```

Block code “fences”

```
Sample text here...
```

Syntax highlighting

```
var foo = function (bar) {  
  return bar++;  
};  
  
console.log(foo(5));
```

14.7 Links

link text

link with title

14.8 Images



Minion



Like links, Images also have a footnote style syntax



Alt

text

With a reference later in the document defining the URL location:

INDICES AND TABLES

- `genindex`
- `modindex`
- `search`

PYTHON MODULE INDEX

V

- `vent.io.devices`, [55](#)
- `vent.io.devices.base`, [31](#)
- `vent.io.devices.pins`, [40](#)
- `vent.io.devices.sensors`, [43](#)
- `vent.io.devices.valves`, [49](#)

Symbols

<code>_CONFIG_FIELDS</code> (<i>vent.io.devices.base.ADS1015 attribute</i>), 39	<code>(vent.io.devices.pins.PWMOutput attribute)</code> , 42
<code>_CONFIG_FIELDS</code> (<i>vent.io.devices.base.ADS1115 attribute</i>), 37	<code>_DEFAULT_STORED_OBSERVATIONS</code> (<i>vent.io.devices.sensors.Sensor attribute</i>), 43
<code>_CONFIG_VALUES</code> (<i>vent.io.devices.base.ADS1015 attribute</i>), 39	<code>_DEFAULT_VALUES</code> (<i>vent.io.devices.base.ADS1015 attribute</i>), 39
<code>_CONFIG_VALUES</code> (<i>vent.io.devices.base.ADS1115 attribute</i>), 37	<code>_DEFAULT_VALUES</code> (<i>vent.io.devices.base.ADS1115 attribute</i>), 37
<code>_ControlModuleBase__analyze_last_waveform()</code> (<i>vent.controller.control_module.ControlModuleBase method</i>), 15	<code>_DEFAULT_offset_voltage</code> (<i>vent.io.devices.sensors.AnalogSensor attribute</i>), 45
<code>_ControlModuleBase__calculate_control_signal_in()</code> (<i>vent.controller.control_module.ControlModuleBase method</i>), 15	<code>_DEFAULT_output_span</code> (<i>vent.io.devices.sensors.AnalogSensor attribute</i>), 45
<code>_ControlModuleBase__get_PID_error()</code> (<i>vent.controller.control_module.ControlModuleBase method</i>), 15	<code>_FLOW_OFFSET</code> (<i>vent.io.devices.sensors.SFM3200 attribute</i>), 48
<code>_ControlModuleBase__test_critical_levels()</code> (<i>vent.controller.control_module.ControlModuleBase method</i>), 15	<code>_FLOW_SCALE_FACTOR</code> (<i>vent.io.devices.sensors.SFM3200 attribute</i>), 48
<code>_ControlModuleBase__update_alarms()</code> (<i>vent.controller.control_module.ControlModuleBase method</i>), 15	<code>_FORMS</code> (<i>vent.io.devices.valves.OnOffValve attribute</i>), 51
<code>_ControlModuleSimulator__SimulatedPropValve()</code> (<i>vent.controller.control_module.ControlModuleSimulator method</i>), 16	<code>_FORMS</code> (<i>vent.io.devices.valves.SolenoidBase attribute</i>), 50
<code>_ControlModuleSimulator__SimulatedSolenoid()</code> (<i>vent.controller.control_module.ControlModuleSimulator method</i>), 16	<code>_HARDWARE_PWM_PINS</code> (<i>vent.io.devices.pins.PWMOutput attribute</i>), 42
<code>_DEFAULT_ADDRESS</code> (<i>vent.io.devices.base.ADS1015 attribute</i>), 39	<code>PID_reset()</code> (<i>vent.controller.control_module.ControlModuleBase method</i>), 15
<code>_DEFAULT_ADDRESS</code> (<i>vent.io.devices.base.ADS1115 attribute</i>), 37	<code>PID_update()</code> (<i>vent.controller.control_module.ControlModuleBase method</i>), 15
<code>_DEFAULT_ADDRESS</code> (<i>vent.io.devices.sensors.SFM3200 attribute</i>), 48	<code>_PIGPIO_MODES</code> (<i>vent.io.devices.pins.Pin attribute</i>), 40
<code>_DEFAULT_CALIBRATION</code> (<i>vent.io.devices.sensors.AnalogSensor attribute</i>), 45	<code>_POINTER_FIELDS</code> (<i>vent.io.devices.base.ADS1015 attribute</i>), 39
<code>_DEFAULT_FREQUENCY</code> (<i>vent.io.devices.pins.PWMOutput attribute</i>), 42	<code>_POINTER_FIELDS</code> (<i>vent.io.devices.base.ADS1115 attribute</i>), 37
<code>_DEFAULT_SOFT_FREQ</code>	<code>_POINTER_VALUES</code> (<i>vent.io.devices.base.ADS1015 attribute</i>), 39
	<code>_POINTER_VALUES</code> (<i>vent.io.devices.base.ADS1115 attribute</i>), 37
	<code>_PWMOutput__hardware_pwm()</code> (<i>vent.io.devices.pins.PWMOutput method</i>), 42

```

_PWMOutput__pwm()
    (vent.io.devices.pins.PWMOutput method), 43
_PWMOutput__software_pwm()
    (vent.io.devices.pins.PWMOutput method), 43
_TIMEOUT (vent.io.devices.base.ADS1115 attribute), 37
__init__() (vent.alarm.alarm.Alarm method), 63
__init__() (vent.alarm.condition.AlarmSeverityCondition
    method), 68
__init__() (vent.alarm.condition.Condition method), 65
__init__() (vent.alarm.condition.TimeValueCondition
    method), 67
__init__() (vent.alarm.condition.ValueCondition
    method), 66
__init__() (vent.common.message.ControlSetting
    method), 8
__init__() (vent.common.message.SensorValues
    method), 8
__init__() (vent.common.values.Value method), 10
__init__() (vent.coordinator.coordinator.CoordinatorLocal
    method), 18
__init__() (vent.gui.widgets.monitor.Monitor
    method), 25
__init__() (vent.gui.widgets.status_bar.HeartBeat
    method), 28
__init__() (vent.io.devices.base.ADS1015 method), 39
__init__() (vent.io.devices.base.ADS1115 method), 37
__init__() (vent.io.devices.base.I2CDevice method), 33
__init__() (vent.io.devices.base.I2CDevice.Register
    method), 34
__init__() (vent.io.devices.base.I2CDevice.Register.Value
    method), 35
__init__() (vent.io.devices.base.IODeviceBase
    method), 32
__init__() (vent.io.devices.base.SPIDevice method), 36
__init__() (vent.io.devices.pins.PWMOutput
    method), 42
__init__() (vent.io.devices.pins.Pin method), 40
__init__() (vent.io.devices.sensors.AnalogSensor
    method), 45
__init__() (vent.io.devices.sensors.SFM3200
    method), 48
__init__() (vent.io.devices.sensors.Sensor method), 44
__init__() (vent.io.devices.sensors.SimSensor
    method), 49
__init__() (vent.io.devices.valves.OnOffValve
    method), 51
__init__() (vent.io.devices.valves.PWMControlValve
    method), 52
__init__() (vent.io.devices.valves.SimControlValve
    method), 54
__init__() (vent.io.devices.valves.SolenoidBase
    method), 50
__init__() (vent.io.hal.Hal method), 56
_abc_impl (vent.io.devices.sensors.AnalogSensor at-
    tribute), 46
_abc_impl (vent.io.devices.sensors.DLiteSensor at-
    tribute), 47
_abc_impl (vent.io.devices.sensors.SFM3200 at-
    tribute), 48
_abc_impl (vent.io.devices.sensors.Sensor attribute), 45
_abc_impl (vent.io.devices.sensors.SimSensor at-
    tribute), 49
_abc_impl (vent.io.devices.valves.OnOffValve at-
    tribute), 51
_abc_impl (vent.io.devices.valves.PWMControlValve
    attribute), 53
_abc_impl (vent.io.devices.valves.SimControlValve at-
    tribute), 55
_abc_impl (vent.io.devices.valves.SimOnOffValve at-
    tribute), 54
_abc_impl (vent.io.devices.valves.SolenoidBase
    attribute), 50
_alarm_to_COPY() (vent.controller.control_module.ControlModuleBase
    method), 14
_check_and_set_attr()
    (vent.io.devices.sensors.AnalogSensor
    method), 46
_child (vent.alarm.condition.Condition attribute), 65
_clear() (vent.io.devices.sensors.Sensor method), 44
_close() (vent.io.devices.base.I2CDevice method), 33
_close() (vent.io.devices.base.IODeviceBase
    method), 32
_close() (vent.io.devices.base.SPIDevice method), 36
_controls_from_COPY()
    (vent.controller.control_module.ControlModuleBase
    method), 14
_convert() (vent.io.devices.sensors.AnalogSensor
    method), 46
_convert() (vent.io.devices.sensors.DLiteSensor
    method), 47
_convert() (vent.io.devices.sensors.SFM3200
    method), 48
_convert() (vent.io.devices.sensors.Sensor method), 44
_convert() (vent.io.devices.sensors.SimSensor
    method), 49
_duty() (vent.io.devices.pins.PWMOutput method), 42
_fill_attr() (vent.io.devices.sensors.AnalogSensor
    method), 46

```

_get_control_signal_in() (vent.controller.control_module.ControlModuleBase method), 14
 _get_control_signal_out() (vent.controller.control_module.ControlModuleBase method), 15
 _heartbeat() (vent.gui.widgets.status_bar.HeartBeat method), 28
 _initialize_set_to_COPY() (vent.controller.control_module.ControlModuleBase method), 14
 _instance(vent.alarm.alarm_manager.Alarm_Manager attribute), 61
 _is_running(vent.coordinator.coordinator.CoordinatorLocal attribute), 18
 _limits_changed() (vent.gui.widgets.monitor.Monitor method), 25
 _load_valve_response() (vent.io.devices.valves.PWMControlValve method), 53
 _mid_check(vent.alarm.condition.CycleAlarmSeverityCondition attribute), 69
 _mid_check(vent.alarm.condition.CycleValueCondition attribute), 66
 _open() (vent.io.devices.base.I2CDevice method), 33
 _open() (vent.io.devices.base.SPIDevice method), 36
 _raw_read() (vent.io.devices.sensors.AnalogSensor method), 46
 _raw_read() (vent.io.devices.sensors.SFM3200 method), 48
 _raw_read() (vent.io.devices.sensors.Sensor method), 44
 _raw_read() (vent.io.devices.sensors.SimSensor method), 49
 _read() (vent.io.devices.sensors.AnalogSensor method), 46
 _read() (vent.io.devices.sensors.Sensor method), 44
 _read_conversion() (vent.io.devices.base.ADS1115 method), 38
 _read_last_cfg() (vent.io.devices.base.ADS1115 method), 39
 _ready() (vent.io.devices.base.ADS1115 method), 39
 _reset() (vent.controller.control_module.Balloon_Simulator method), 16
 _safe_limits_changed() (vent.gui.widgets.plot.Plot method), 26
 _sensor_to_COPY() (vent.controller.control_module.ControlModuleBase method), 14
 _sensor_to_COPY() (vent.controller.control_module.ControlModuleDevice method), 15
 _sensor_to_COPY() (vent.controller.control_module.ControlModuleSimulator method), 16
 _start() (vent.io.devices.sensors.SFM3200 method), 48
 _start_cycle(vent.alarm.condition.CycleAlarmSeverityCondition attribute), 68
 _start_cycle(vent.alarm.condition.CycleValueCondition attribute), 66
 _start_mainloop() (vent.controller.control_module.ControlModuleBase method), 15
 _start_mainloop() (vent.controller.control_module.ControlModuleDevice method), 16
 _start_mainloop() (vent.controller.control_module.ControlModuleSimulator method), 16
 _verify() (vent.io.devices.sensors.AnalogSensor method), 46
 _verify() (vent.io.devices.sensors.SFM3200 method), 48
 _verify() (vent.io.devices.sensors.Sensor method), 44
 _verify() (vent.io.devices.sensors.SimSensor method), 49

A

abs_range() (vent.common.values.Value property), 10
 active_alarms (in module vent.alarm.alarm_manager), 60
 active_alarms(vent.alarm.alarm_manager.Alarm_Manager attribute), 60, 61
 add_callback() (vent.alarm.alarm_manager.Alarm_Manager method), 62
 additional_values (vent.common.message.SensorValues attribute), 8
 ADS1015 (class in vent.io.devices.base), 39
 ADS1115 (class in vent.io.devices.base), 36
 age() (vent.io.devices.sensors.Sensor method), 44
 Alarm (class in vent.alarm.alarm), 63
 alarm(vent.gui.widgets.monitor.Monitor attribute), 25
 alarm_level() (vent.gui.widgets.status_bar.Message_Display property), 27
 Alarm_Manager (class in vent.alarm.alarm_manager), 60
 Alarm_Rule (class in vent.alarm.rule), 69
 alarm_state() (vent.gui.widgets.monitor.Monitor property), 25
 alarm_type() (vent.alarm.alarm.Alarm property), 64
 AlarmManager (class in vent.gui.alarm_manager), 29
 AlarmSeverity (class in vent.alarm), 59
 AlarmSeverityCondition (class in vent.alarm.condition), 67

AlarmType (class in *vent.alarm*), 59

AnalogSensor (class in *vent.io.devices.sensors*), 45

aux_pressure() (*vent.io.hal.Hal* property), 57

B

Balloon_Simulator (class in *vent.controller.control_module*), 16

bel6_to_native() (in module *vent.io.devices.base*), 39

beatheart() (*vent.gui.widgets.status_bar.HeartBeat* method), 28

BREATHS_PER_MINUTE (*vent.common.values.ValueName* attribute), 9

C

calibrate() (*vent.io.devices.sensors.AnalogSensor* method), 46

calibrate() (*vent.io.devices.sensors.DLiteSensor* method), 47

callbacks (*vent.alarm.alarm_manager.Alarm_Manager* attribute), 61

cfg() (*vent.io.devices.base.ADS1115* property), 38

check() (*vent.alarm.condition.AlarmSeverityCondition* method), 68

check() (*vent.alarm.condition.Condition* method), 65

check() (*vent.alarm.condition.CycleAlarmSeverityCondition* method), 69

check() (*vent.alarm.condition.CycleValueCondition* method), 67

check() (*vent.alarm.condition.TimeValueCondition* method), 67

check() (*vent.alarm.condition.ValueCondition* method), 66

check() (*vent.alarm.rule.Alarm_Rule* method), 70

check_alarm() (*vent.gui.widgets.monitor.Monitor* method), 25

check_rule() (*vent.alarm.alarm_manager.Alarm_Manager* method), 61

check_timeout() (*vent.gui.widgets.status_bar.HeartBeat* method), 28

clear_all_alarms() (*vent.alarm.alarm_manager.Alarm_Manager* method), 62

clear_message() (*vent.gui.widgets.status_bar.Message_Display* method), 27

cleared_alarms (*vent.alarm.alarm_manager.Alarm_Manager* attribute), 61

close() (*vent.io.devices.valves.OnOffValve* method), 51

close() (*vent.io.devices.valves.PWMControlValve* method), 52

close() (*vent.io.devices.valves.SimControlValve* method), 55

close() (*vent.io.devices.valves.SimOnOffValve* method), 54

close() (*vent.io.devices.valves.SolenoidBase* method), 50

Condition (class in *vent.alarm.condition*), 64

config() (*vent.io.devices.base.ADS1115* property), 38

CONTROL (in module *vent.common.values*), 11

control() (*vent.common.values.Value* property), 10

ControlModuleBase (class in *vent.controller.control_module*), 13

ControlModuleDevice (class in *vent.controller.control_module*), 15

ControlModuleSimulator (class in *vent.controller.control_module*), 16

ControlSetting (class in *vent.common.message*), 8

CoordinatorBase (class in *vent.coordinator.coordinator*), 17

CoordinatorLocal (class in *vent.coordinator.coordinator*), 17

CoordinatorRemote (class in *vent.coordinator.coordinator*), 18

CycleAlarmSeverityCondition (class in *vent.alarm.condition*), 68

CycleValueCondition (class in *vent.alarm.condition*), 66

D

data() (*vent.io.devices.sensors.Sensor* property), 44

deactivate() (*vent.alarm.alarm.Alarm* method), 64

deactivate_alarm() (*vent.alarm.alarm_manager.Alarm_Manager* method), 62

decimals() (*vent.common.values.Value* property), 10

default() (*vent.common.values.Value* property), 10

dependencies (*vent.alarm.alarm_manager.Alarm_Manager* attribute), 61

dismiss_alarm() (*vent.alarm.alarm_manager.Alarm_Manager* method), 62

DLiteSensor (class in *vent.io.devices.sensors*), 46

do_pid_control() (*vent.controller.control_module.ControlModuleBase* method), 15

do_state_control() (*vent.controller.control_module.ControlModuleBase* method), 15

draw_state() (*vent.gui.widgets.status_bar.Message_Display* method), 27

duty() (*vent.io.devices.pins.PWMOutput* property), 42

E

emit_alarm() (*vent.alarm.alarm_manager.Alarm_Manager* method), 61

Error (class in *vent.common.message*), 8

extract() (*vent.io.devices.base.I2CDevice.Register.ValueField* method), 35

F

FIO2 (*vent.common.values.ValueName* attribute), 9
 flow_ex() (*vent.io.hal.Hal* property), 57
 flow_in() (*vent.io.hal.Hal* property), 57
 form() (*vent.io.devices.valves.SolenoidBase* property), 50
 frequency() (*vent.io.devices.pins.PWMOutput* property), 42

G

get() (*vent.io.devices.sensors.Sensor* method), 44
 get_active_alarms() (*vent.controller.control_module.ControlModuleBase* method), 14
 get_alarm_manager() (in module *vent.alarm.condition*), 64
 get_alarm_manager() (in module *vent.gui.alarm_manager*), 29
 get_alarm_severity() (*vent.alarm.alarm_manager.Alarm_Manager* method), 62
 get_control() (*vent.controller.control_module.ControlModuleBase* method), 14
 get_control() (*vent.coordinator.coordinator.CoordinatorBase* method), 17
 get_control() (*vent.coordinator.coordinator.CoordinatorLocal* method), 18
 get_control() (*vent.coordinator.coordinator.CoordinatorRemote* method), 18
 get_control_module() (in module *vent.controller.control_module*), 16
 get_coordinator() (in module *vent.coordinator.coordinator*), 19
 get_logged_alarms() (*vent.controller.control_module.ControlModuleBase* method), 14
 get_past_waveforms() (*vent.controller.control_module.ControlModuleBase* method), 15
 get_pressure() (*vent.controller.control_module.Balloon_Simulation* method), 16
 get_sensors() (*vent.controller.control_module.ControlModuleBase* method), 14
 get_sensors() (*vent.coordinator.coordinator.CoordinatorBase* method), 17
 get_sensors() (*vent.coordinator.coordinator.CoordinatorDefault* method), 18
 get_sensors() (*vent.coordinator.coordinator.CoordinatorRemote* method), 18
 get_target_waveform() (*vent.controller.control_module.ControlModuleBase* method), 15
 get_volume() (*vent.controller.control_module.Balloon_Simulation* method), 16

H

Hal (class in *vent.io.hal*), 55
 handle() (*vent.io.devices.base.IODeviceBase* property), 32
 hardware_enabled() (*vent.io.devices.pins.PWMOutput* property), 42
 HeartBeat (class in *vent.gui.widgets.status_bar*), 27
 heartbeat (*vent.gui.widgets.status_bar.HeartBeat* attribute), 28
 heartbeat() (*vent.coordinator.process_manager.ProcessManager* method), 19
 HIGH (*vent.alarm.AlarmSeverity* attribute), 59
 HIGH_O2 (*vent.alarm.AlarmType* attribute), 59
 HIGH_PEEP (*vent.alarm.AlarmType* attribute), 59
 HIGH_PRESSURE (*vent.alarm.AlarmType* attribute), 59
 HIGH_VTE (*vent.alarm.AlarmType* attribute), 59
 HUMIDITY (*vent.common.values.ValueName* attribute), 9

I

I2CDevice (class in *vent.io.devices.base*), 32
 I2CDevice.Register (class in *vent.io.devices.base*), 34
 I2CDevice.Register.ValueField (class in *vent.io.devices.base*), 34
 Incrementer (*vent.alarm.alarm.Alarm* attribute), 63
 IE_RATIO (*vent.common.values.ValueName* attribute), 9
 init_ui() (*vent.gui.widgets.monitor.Monitor* method), 25
 init_ui() (*vent.gui.widgets.status_bar.HeartBeat* method), 28
 init_ui() (*vent.gui.widgets.status_bar.Message_Display* method), 27
 init_ui() (*vent.gui.widgets.status_bar.Power_Button* method), 28
 init_ui() (*vent.gui.widgets.status_bar.Status_Bar* method), 26
 insert() (*vent.io.devices.base.I2CDevice.Register.ValueField* method), 35
 INSPIRATION_TIME_SEC (*vent.common.values.ValueName* attribute), 9
 inverse_response() (*vent.io.devices.valves.PWMControlValve* method), 53
 IODeviceBase (class in *vent.io.devices.base*), 31
 is_open() (*vent.io.devices.valves.OnOffValve* property), 51
 is_open() (*vent.io.devices.valves.PWMControlValve* property), 52
 is_open() (*vent.io.devices.valves.SimControlValve* property), 54
 is_open() (*vent.io.devices.valves.SimOnOffValve* property), 54

- [is_open\(\)](#) (*vent.io.devices.valves.SolenoidBase* property), 50
[is_running\(\)](#) (*vent.controller.control_module.ControlModuleBase* method), 15
[is_running\(\)](#) (*vent.coordinator.coordinator.CoordinatorBase* method), 17
[is_running\(\)](#) (*vent.coordinator.coordinator.CoordinatorLocal* method), 18
[is_running\(\)](#) (*vent.coordinator.coordinator.CoordinatorRemote* method), 19
- ## L
- [LEAK](#) (*vent.alarm.AlarmType* attribute), 59
[level_changed](#) (*vent.gui.widgets.status_bar.Message_Display* attribute), 27
[LIMITS](#) (in module *vent.common.values*), 11
[limits_changed](#) (*vent.gui.widgets.monitor.Monitor* attribute), 25
[limits_changed](#) (*vent.gui.widgets.plot.Plot* attribute), 26
[load_rule\(\)](#) (*vent.alarm.alarm_manager.Alarm_Manager* method), 61
[load_rules\(\)](#) (*vent.alarm.alarm_manager.Alarm_Manager* method), 61
[logged_alarms](#) (*vent.alarm.alarm_manager.Alarm_Manager* attribute), 61
[LOW](#) (*vent.alarm.AlarmSeverity* attribute), 60
[LOW_O2](#) (*vent.alarm.AlarmType* attribute), 59
[LOW_PEEP](#) (*vent.alarm.AlarmType* attribute), 59
[LOW_PRESSURE](#) (*vent.alarm.AlarmType* attribute), 59
[LOW_VTE](#) (*vent.alarm.AlarmType* attribute), 59
- ## M
- [make_icons\(\)](#) (*vent.gui.widgets.status_bar.Message_Display* method), 27
[manager](#) (*vent.alarm.condition.Condition* attribute), 64
[manager\(\)](#) (*vent.alarm.alarm.Alarm* property), 64
[manager\(\)](#) (*vent.alarm.condition.Condition* property), 65
[maxlen_data\(\)](#) (*vent.io.devices.sensors.Sensor* property), 44
[MEDIUM](#) (*vent.alarm.AlarmSeverity* attribute), 60
[message_cleared](#) (*vent.gui.widgets.status_bar.Message_Display* attribute), 27
[Message_Display](#) (class in *vent.gui.widgets.status_bar*), 26
[mode\(\)](#) (*vent.alarm.condition.AlarmSeverityCondition* property), 68
[mode\(\)](#) (*vent.alarm.condition.ValueCondition* property), 66
[mode\(\)](#) (*vent.io.devices.pins.Pin* property), 41
[module](#)
 vent.alarm, 59
 vent.alarm.alarm, 63
 vent.alarm.alarm_manager, 60
 vent.alarm.condition, 64
 vent.alarm.rule, 69
 vent.common.message, 7
 vent.common.values, 8
 vent.controller.control_module, 13
 vent.coordinator.coordinator, 17
 vent.coordinator.process_manager, 19
 vent.gui.alarm_manager, 29
 vent.gui.styles, 28
 vent.gui.widgets.components, 28
 vent.gui.widgets.control, 24
 vent.gui.widgets.monitor, 24
 vent.gui.widgets.plot, 25
 vent.gui.widgets.status_bar, 26
 vent.io, 57
 vent.io.devices, 55
 vent.io.devices.base, 31
 vent.io.devices.pins, 40
 vent.io.devices.sensors, 43
 vent.io.devices.valves, 49
 vent.io.hal, 55
 Monitor (class in *vent.gui.widgets.monitor*), 24
 monitor_alarm() (*vent.gui.alarm_manager.AlarmManager* method), 29
 MONITOR_UPDATE_INTERVAL (in module *vent.gui.styles*), 28
- ## N
- [n_cycles\(\)](#) (*vent.alarm.condition.CycleAlarmSeverityCondition* property), 69
[n_cycles\(\)](#) (*vent.alarm.condition.CycleValueCondition* property), 66
[name\(\)](#) (*vent.common.values.Value* property), 10
[native16_to_be\(\)](#) (in module *vent.io.devices.base*), 40
[new_alarm](#) (*vent.gui.alarm_manager.AlarmManager* attribute), 29
- ## O
- [OBSTRUCTION](#) (*vent.alarm.AlarmType* attribute), 59
[OFF](#) (*vent.alarm.AlarmSeverity* attribute), 60
[OnOffValve](#) (class in *vent.io.devices.valves*), 50
[open\(\)](#) (*vent.io.devices.valves.OnOffValve* method), 51
[open\(\)](#) (*vent.io.devices.valves.PWMControlValve* method), 52
[open\(\)](#) (*vent.io.devices.valves.SimControlValve* method), 55
[open\(\)](#) (*vent.io.devices.valves.SimOnOffValve* method), 53
[open\(\)](#) (*vent.io.devices.valves.SolenoidBase* method), 50
[OUpdate\(\)](#) (*vent.controller.control_module.Balloon_Simulator* method), 16

P

[pack \(\) \(vent.io.devices.base.I2CDevice.Register method\), 34](#)
[pack \(\) \(vent.io.devices.base.I2CDevice.Register.ValueField method\), 35](#)
[parse_message \(\) \(vent.gui.alarm_manager.AlarmManager method\), 29](#)
[PEEP \(vent.common.values.ValueName attribute\), 9](#)
[PEEP_TIME \(vent.common.values.ValueName attribute\), 9](#)
[pending_clears \(vent.alarm.alarm_manager.AlarmManager attribute\), 61](#)
[pig \(\) \(vent.io.devices.base.IODeviceBase property\), 32](#)
[pigpiod_ok \(\) \(vent.io.devices.base.IODeviceBase property\), 32](#)
[Pin \(class in vent.io.devices.pins\), 40](#)
[PIP \(vent.common.values.ValueName attribute\), 8](#)
[PIP_TIME \(vent.common.values.ValueName attribute\), 9](#)
[Plot \(class in vent.gui.widgets.plot\), 25](#)
[PLOT_FREQ \(in module vent.gui.widgets.plot\), 25](#)
[PLOT_TIMER \(in module vent.gui.widgets.plot\), 25](#)
[Power_Button \(class in vent.gui.widgets.status_bar\), 28](#)
[PRESSURE \(vent.common.values.ValueName attribute\), 9](#)
[pressure \(\) \(vent.io.hal.Hal property\), 57](#)
[print_config \(\) \(vent.io.devices.base.ADS1115 method\), 38](#)
[ProcessManager \(class in vent.coordinator.process_manager\), 19](#)
[PWMControlValve \(class in vent.io.devices.valves\), 51](#)
[PWMOutput \(class in vent.io.devices.pins\), 41](#)

R

[read \(\) \(vent.io.devices.pins.Pin method\), 41](#)
[read \(\) \(vent.io.devices.pins.PWMOutput method\), 42](#)
[read_conversion \(\) \(vent.io.devices.base.ADS1115 method\), 38](#)
[read_device \(\) \(vent.io.devices.base.I2CDevice method\), 33](#)
[read_register \(\) \(vent.io.devices.base.I2CDevice method\), 33](#)
[register_alarm \(\) \(vent.alarm.alarm_manager.AlarmManager method\), 28](#)
[register_dependency \(\) \(vent.alarm.alarm_manager.AlarmManager method\), 62](#)
[reset \(\) \(vent.alarm.alarm_manager.AlarmManager method\), 63](#)
[reset \(\) \(vent.alarm.condition.AlarmSeverityCondition method\), 68](#)
[reset \(\) \(vent.alarm.condition.Condition method\), 65](#)
[reset \(\) \(vent.alarm.condition.CycleAlarmSeverityCondition method\), 69](#)
[reset \(\) \(vent.alarm.condition.CycleValueCondition method\), 67](#)
[reset \(\) \(vent.alarm.condition.TimeValueCondition method\), 67](#)
[reset \(\) \(vent.alarm.condition.ValueCondition method\), 66](#)
[reset \(\) \(vent.alarm.rule.Alarm_Rule method\), 70](#)
[reset \(\) \(vent.io.devices.sensors.Sensor method\), 44](#)
[reset \(\) \(vent.io.devices.sensors.SFM3200 method\), 48](#)
[response \(\) \(vent.io.devices.valves.PWMControlValve method\), 52](#)
[restart_process \(\) \(vent.coordinator.process_manager.ProcessManager method\), 19](#)
[rules \(vent.alarm.alarm_manager.Alarm_Manager attribute\), 61](#)

S

[safe_range \(\) \(vent.common.values.Value property\), 10](#)
[Sensor \(class in vent.io.devices.sensors\), 43](#)
[SENSOR \(in module vent.common.values\), 10](#)
[sensor \(\) \(vent.common.values.Value property\), 10](#)
[SensorValues \(class in vent.common.message\), 7](#)
[set_alarm \(\) \(vent.gui.widgets.monitor.Monitor method\), 25](#)
[set_control \(\) \(vent.controller.control_module.ControlModuleBase method\), 14](#)
[set_control \(\) \(vent.coordinator.coordinator.CoordinatorBase method\), 17](#)
[set_control \(\) \(vent.coordinator.coordinator.CoordinatorLocal method\), 18](#)
[set_control \(\) \(vent.coordinator.coordinator.CoordinatorRemote method\), 18](#)
[set_dark_palette \(\) \(in module vent.gui.styles\), 29](#)
[set_duration \(\) \(vent.gui.widgets.plot.Plot method\), 26](#)
[set_flow_in \(\) \(vent.controller.control_module.Balloon_Simulator method\), 16](#)
[set_flow_out \(\) \(vent.controller.control_module.Balloon_Simulator method\), 16](#)
[set_indicator \(\) \(vent.gui.widgets.status_bar.HeartBeat method\), 28](#)
[set_safe_limits \(\) \(vent.gui.widgets.plot.Plot method\), 26](#)
[setpoint \(\) \(vent.io.devices.valves.PWMControlValve property\), 52](#)
[setpoint \(\) \(vent.io.devices.valves.SimControlValve property\), 55](#)
[setpoint_ex \(\) \(vent.io.hal.Hal property\), 57](#)
[setpoint_in \(\) \(vent.io.hal.Hal property\), 57](#)
[severity \(\) \(vent.alarm.alarm.Alarm property\), 64](#)

severity() (*vent.alarm.rule.Alarm_Rule* property), 70
 SFM3200 (*class in vent.io.devices.sensors*), 47
 SimControlValve (*class in vent.io.devices.valves*), 54
 SimOnOffValve (*class in vent.io.devices.valves*), 53
 SimSensor (*class in vent.io.devices.sensors*), 48
 snoozed_alarms (*vent.alarm.alarm_manager.Alarm_Manager* attribute), 61
 SolenoidBase (*class in vent.io.devices.valves*), 49
 SPIDevice (*class in vent.io.devices.base*), 35
 start() (*vent.controller.control_module.ControlModuleBase* method), 15
 start() (*vent.coordinator.coordinator.CoordinatorBase* method), 17
 start() (*vent.coordinator.coordinator.CoordinatorLocal* method), 18
 start() (*vent.coordinator.coordinator.CoordinatorRemote* method), 18
 start_process() (*vent.coordinator.process_manager.ProcessManager* method), 19
 start_timer() (*vent.gui.widgets.status_bar.HeartBeat* method), 28
 staticMetaObject (*vent.gui.alarm_manager.AlarmManager* attribute), 29
 staticMetaObject (*vent.gui.widgets.monitor.Monitor* attribute), 25
 staticMetaObject (*vent.gui.widgets.plot.Plot* attribute), 26
 staticMetaObject (*vent.gui.widgets.status_bar.HeartBeat* attribute), 28
 staticMetaObject (*vent.gui.widgets.status_bar.Message_Display* attribute), 27
 staticMetaObject (*vent.gui.widgets.status_bar.Power_Button* attribute), 28
 staticMetaObject (*vent.gui.widgets.status_bar.Status_Bar* attribute), 26
 Status_Bar (*class in vent.gui.widgets.status_bar*), 26
 stop() (*vent.controller.control_module.ControlModuleBase* method), 15
 stop() (*vent.coordinator.coordinator.CoordinatorBase* method), 17
 stop() (*vent.coordinator.coordinator.CoordinatorLocal* method), 18
 stop() (*vent.coordinator.coordinator.CoordinatorRemote* method), 19
 stop_timer() (*vent.gui.widgets.status_bar.HeartBeat* method), 28
 timeout (*vent.gui.widgets.status_bar.HeartBeat* attribute), 28
 TimeValueCondition (*class in vent.alarm.condition*), 67
 to_dict() (*vent.common.message.SensorValues* method), 8
 to_dict() (*vent.common.values.Value* method), 10
 toggle() (*vent.io.devices.pins.Pin* method), 41
 toggle_alarm() (*vent.gui.widgets.monitor.Monitor* method), 25
 toggle_control() (*vent.gui.widgets.monitor.Monitor* method), 25
 try_stop_process() (*vent.coordinator.process_manager.ProcessManager* method), 19

U

unpack() (*vent.io.devices.base.I2CDevice.Register* method), 34
 unpack() (*vent.io.devices.base.I2CDevice.Register.ValueField* method), 35
 update() (*vent.alarm.alarm_manager.Alarm_Manager* method), 61
 update() (*vent.controller.control_module.Balloon_Simulator* method), 16
 update() (*vent.io.devices.sensors.Sensor* method), 44
 update_alarms() (*vent.gui.alarm_manager.AlarmManager* method), 29
 update_boxes() (*vent.gui.widgets.monitor.Monitor* method), 25
 update_dependencies() (*vent.alarm.alarm_manager.Alarm_Manager* method), 62
 update_limits() (*vent.gui.widgets.monitor.Monitor* method), 25
 update_message() (*vent.gui.widgets.status_bar.Message_Display* method), 27
 update_value() (*vent.gui.widgets.monitor.Monitor* method), 25
 update_value() (*vent.gui.widgets.plot.Plot* method), 26
 USER_CONFIGURABLE_FIELDS (*vent.io.devices.base.ADS1015* attribute), 39
 USER_CONFIGURABLE_FIELDS (*vent.io.devices.base.ADS1115* attribute), 37

V

Value (*class in vent.common.values*), 9
 ValueCondition (*class in vent.alarm.condition*), 65
 ValueName (*class in vent.common.values*), 8
 vent.alarm module, 59

T

TEMP (*vent.common.values.ValueName* attribute), 9
 timed_update() (*vent.gui.widgets.monitor.Monitor* method), 25

vent.alarm.alarm
 module, 63
 vent.alarm.alarm_manager
 module, 60
 vent.alarm.condition
 module, 64
 vent.alarm.rule
 module, 69
 vent.common.message
 module, 7
 vent.common.values
 module, 8
 vent.controller.control_module
 module, 13
 vent.coordinator.coordinator
 module, 17
 vent.coordinator.process_manager
 module, 19
 vent.gui.alarm_manager
 module, 29
 vent.gui.styles
 module, 28
 vent.gui.widgets.components
 module, 28
 vent.gui.widgets.control
 module, 24
 vent.gui.widgets.monitor
 module, 24
 vent.gui.widgets.plot
 module, 25
 vent.gui.widgets.status_bar
 module, 26
 vent.io
 module, 57
 vent.io.devices
 module, 55
 vent.io.devices.base
 module, 31
 vent.io.devices.pins
 module, 40
 vent.io.devices.sensors
 module, 43
 vent.io.devices.valves
 module, 49
 vent.io.hal
 module, 55
 VTE (*vent.common.values.ValueName attribute*), 9

W

write() (*vent.io.devices.pins.Pin method*), 41
 write() (*vent.io.devices.pins.PWMOutput method*), 42
 write_device() (*vent.io.devices.base.I2CDevice method*), 33